


# Pharma-Jelly Net



## Description:

Pharma-Jelly Net is a Lipo-Colloid non-adherent wound dressings. Pharma-Jelly Net is composed of polyester mesh impregnated with superabsorbent powder and petrolatum.

## Features:

•Healing in moist environment:

Hydrocolloid particles gel upon absorption of wound exudates, creating moist environment in favour of wound healing.

• Painless and atraumatic removal:

Due to its fatty, yet non-greasy, nature, together with the gel formation property upon exudates contact, Lipocolloid dressings' removal is painless and atraumatic for the wound.

• Promote fast healing by stimulating fibroblasts proliferation.

• Highly conformable:

The dressing is very soft and conformable, so it conforms to wound contours.

## Intended Purpose:

Pharma-Jelly Net is indicated for management of non to low exuding wounds.

It could be used for: burns, traumatic wounds, abrasions, post-operative wounds, leg ulcers, pressure ulcers, diabetic foot ulcers, and epidermolysis bullosa.

# Pharma-Jelly Net

Code	Size	Pieces/box	Pieces/Case
JN5050	5 x 5 cm	10	48
JN50100	5 x 10 cm	10	48
		100	16
JN100100	10 x 10 cm	10	48
		100	16
JN100150	10 x 15 cm	10	16
		50	16
JN100200	10 x 20 cm	10	24
		50	16
JN100300	10 x 30 cm	10	24
		30	16
JN150150	15 x 15 cm	10	24
		40	24
JN150200	15 x 20 cm	10	24
		25	24
JN200200	20 x 20 cm	10	24
		25	24
JN200300	20 x 30 cm	10	24
		25	24
JN120250	12 x 25 cm	10	24
		25	24
JN150250	15 x 25 cm	10	24
		25	24
JN150300	15 x 30 cm	10	24
		25	24
JN7575	7.5 x 7.5 cm	10	48
JN75100	7.5 x 10 cm	10	48

# Pharma-Jelly Net

Code	Size	Pieces/box	Pieces/Case
JN75200	7.5 x 20 cm	10	24
JN75400	7.5 x 40 cm	10	24
		25	24
JN751300	7.5 x 1.3 m	1	48
JN751000	7.5 cm x 1m	1	48
JN1001000	10 cm x 1m	1	48
JN1501000	15 cm x 1m	1	48
JN2001000	20 cm x 1 m	1	48
JN1002000	10 cm x 2 m	1	48
JN1502000	15 cm x 2 m	1	48
JN2002000	20 cm x 2 m	1	48
JN100400	10x40 cm	10	24
JN200400	20x40 cm	10	24
JN501000	5x100 cm	2	48
JN757000	7.5x700 cm	1	48
JN1007000	10x700 cm	1	48
JN1507000	15x700 cm	1	48
JN75150	7.5x15 cm	10	24
JN150400	15x40 cm	10	24

- Amria - Alexandria - Egypt  
 Tel. : 002 03 4500264 Fax : 002 03 4500263  
 - Kafr-El-Zayat P.O.Box 38 Egypt  
 Tel. : 002 040 2545643 Fax : 002 040 2545645  
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**  
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3 )						
Customer : Pharmaplast			Product name : Pharma-Jelly net					
Issue date : xx-xx-xxxx			Prepared by : Marina Medhat					
Colors : CMYK			Dimension : A4					
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	09-02-2022	Updating codes	File maintenance					