

AbsoClear® I.V Dressings CHG



Description:

AbsoClear® I.V dressings are a vast range of products that use the **AbsoClear®** gel technology with chlorhexidine to secure catheters, hoses and pipes. **AbsoClear®** gel is a water clear hydrogel, made of ionic polymer, chlorhexidine and glycerol it is placed as an island pad on a PU film backing coated with hypoallergenic adhesive.

AbsoClear® provides soft contact layer with high absorption capacity.

Features:

- Highly absorbent gel layer.
- Invitro Studies show bacteriostatic and fungistatic effect up to 7 days.
- Tackiness: adheres to skin, so it can be applied without fixation tapes of secondary film dressings.
- Easy removal: it can be removed without pain, trauma or residues on skin.

Benefits:

1. **AbsoClear® I.V** dressings are comprehensively designed to deliver exceptional patient care through advanced catheter securement, superior wear time, gentle removal, and excellent absorption.
2. **AbsoClear® I.V** dressings help secure and preserve the integrity of devices, minimize movement, prevent catheter dislodgement.
3. Invitro tests have shown that **AbsoClear®** exhibits bacteriostatic and fungistatic effect.
4. **AbsoClear® I.V** dressings with CHG are ideal for use during catheter application for prevention of intra-vascular catheter-associated infections

Indications:

- **AbsoClear IV dressings:**

For Cannula and Catheters fixation.

AbsoClear® I.V Dressings CHG

Contraindications:






Any known sensitivity for chlorhexidine

How Supplied:

• **AbsoClear® I.V** dressings are supplied as sterile, individually wrapped in medical sealed paper.


AbsoClear® I.V dressings are sterilized by EO gas.

Sizes

AbsoClear PU IV CHG	Size	Pieces/ box	Boxes / Case	Case dimensions
Abso CHG 76	7 cm x 6 cm	25	24	40 x 32 x 28 cm
Abso CHG 77	7 cm x 7 cm	25	24	40 x 32 x 28 cm
Abso CHG 87	8 cm x 7 cm	25	24	40 x 32 x 28 cm
Abso CHG 810	8 cm x 10 cm	25	24	40 x 32 x 28 cm
Abso CHG 2810	8.25 cm x 10 cm	25	24	40 x 32 x 28 cm
Abso CHG 856	8.5 cm x 6 cm	25	24	40 x 32 x 28 cm
Abso CHG 857	8.5 cm x 7 cm	25	24	40 x 32 x 28 cm
Abso CHG 8510	8.5 cm x 10 cm	25	24	40 x 32 x 28 cm
Abso CHG 97	9 cm x 7 cm	25	24	40 x 32 x 28 cm
Abso CHG 910	9 cm x 10 cm	25	24	40 x 32 x 28 cm
Abso CHG 1012	10 cm x 12 cm	25	24	40 x 32 x 28 cm
Abso CHG 1014	10 m x 14 cm	25	24	40 x 32 x 28 cm
Abso CHG 68	6 cm x 8 cm	25	24	40 x 32 x 28 cm
AbsoClear PU IV FS CHG	Size	Pieces/ box	Boxes / Case	Case dimensions
Abso CHG 657 	6.5 cm x 7 cm	25	24	40 x 32 x 28 cm
Abso CHG 785 	7 cm x 8.5 cm	25	24	40 x 32 x 28 cm
Abso CHG 85115 	8.5 cm x 11.5 cm	25	24	40 x 32 x 28 cm
Abso CHG 1012 	10 cm x 12 cm	25	24	40 x 32 x 28 cm
Abso CHG 10155 	10 cm x 15.5 cm	25	24	40 x 32 x 28 cm

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
www.pharmaplast-online.com

 Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6				
Customer : Pharmaplast			Product name : AbsoClear® I.V Dressings CHG			
Issue date : 25-04-2010			Prepared by : Nader Nabil Zaki			
Colors : CMYK			Dimension : A4			
Issue Number	Revision Date	Description of Change	Rationale For Change	Approved by		
				Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A			