

Cavicol[®]

“Hydroactive Colloid Paste”



Description:

Cavicol[®] is a gelatin free hydrocolloid paste sterilized by radiation. It forms gel upon absorbing wound exudates, providing a moist wound environment and supporting granulation and epithelialization.

Mode of action:

- Paraffinic components have non-polar fatty nature, so they do not adhere to the polar proteinaceous exudates.
- Hydrocolloid particles absorb exudates and turn to gel which provides optimum moist environment that enhances wound healing. Furthermore, it aids non-adherence due to synergistic action between gel and paraffinic components.

Features:

- Healing in moist environment: interacts with wound exudate to form a moist wound environment.
- Painless and atraumatic removal: Due to its fatty, yet non-greasy, nature, together with the gel formation property upon exudates contact, the hydrocolloid paste allows the atraumatic removal of the dressing without damaging newly formed tissue.
- Absorbent: Due to the fact that **Cavicol[®]** contains a high percentage of superabsorbent polymer, the paste aids in the absorption of excessive exudate and extend the life of the secondary dressing, increasing the wear time and decreasing dressing change rate.
- Aids autolytic debridement.
- Helps increase ostomy system wear time and protects the skin.

Intended Purpose:

Cavicol[®] is intended to be used in association with a secondary dressing such as **Protectfilm[®]** for the management of exuding dermal ulcers, including full-thickness wounds, such as leg ulcers, pressure ulcers and diabetic ulcers.

Cavicol[®] also can be used to seal ostomy bags where paste fills deeper cavities and folds, creating a tight seal between the base plate and the skin.

Cavicol [®]	Size	Pieces/ box	Boxes / Case	Case Dimensions
COLLPASTE15	15gm	1	48	
COLLPASTE30	30gm	1	48	
COLLPASTE75	75gm	1	24	
COLLPASTE40	40gm	1	24	
COLLPASTE50	50gm	1	24	
COLLPASTE100	100gm	1	24	
COLLPASTE250	250gm	1	8	

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com


Pharmaplast
www.pharmaplast-online.com

Pharmaplast		Amria 23512, Alexandria - Egypt www.pharmaplast-online.com							Brochure artwork template no: I-TM /AR/001-6 (Version 3)		
Customer : Pharmaplast		Product name : Cavicol									
Issue date : 15-10-2019		Prepared by : Nader Nabil Zaki									
Colors : CMYK		Dimension : A4									
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by					
						Artwork manager	Marketing department	RA Department			
1	N/A	New Issue	N/A								
2	09-02-2022	Updating Codes	File maintenance								
3	24-05-2022	updating table of sizes	File maintenance								