

**First Aid  
Bandages**

**Pharmaplast**



**Heal your**  
Minor Wounds **Fast**  
With Our **Plasters**

- **Cure - Aid® Display**
- **Detectable Plasters**
- **Value Packs**

# First Aid Bandages

Adhesive plasters that consist of an absorbent low-adherent pad placed on different backing materials coated with pressure sensitive adhesive to meet all patient's needs.


## Intended For

covering skin scratches and minor wounds.

## Available with different backing materials



## Available with different pad types




Available with different shapes




Available with different application systems


Wave-Cut application system



Easy open application system



Frame style application system



# Cure-Aid® Display

- Available with different backing materials
- Display of 48 boxes that can be customized as your need
- The range of this pharmacy display can be created under private label design based on minimum order quantity.



## kids Plasters

**Plasters to distract and delight your one with cute different attractive designs**

Cure - Aid® Kids Plasters with low-adherent absorbent pad placed on washproof printed plastic microperforated Provides protection against dirt and physical contaminants.





Product code	Size	20 Pieces/Box
CUREDISKID20	Shapes	10 Characters
	19X63	5 Animals
	19X38	5 Animals

## Cure - Aid® Waterproof

**100 % waterproof plasters suitable for washing, showering, bathing and swimming.**

Cure Aid® Waterproof Plasters with low-adherent absorbent pad placed on PU waterproof, viral proof & bacterial proof material that protect the wound from external contamination.

Product code	Size	20 Pieces/Box
CUREDIS-WTP20	38X38 mm	12
	30X55 mm	8






WATERPROOF



## Cure - Aid® Washproof

**Ideal for working in wet or moist conditions to protect the wound from water.**

Cure - Aid® Washproof Plasters with low-adherent absorbent pad placed on skin color plastic microperforated washproof material.




Product code	Size	20 Pieces/Box
CUREDIS-WP20	19X72 mm	5
	25X72 mm	5
	25Ø Spots	5
	9X38 mm	5

## Cure - Aid® Sheer

**Ideal for working in wet or moist conditions to protect the wound from water.**

Cure - Aid® Sheer Plasters with low-adherent absorbent pad placed on translucent plastic microperforated washproof material to reflect skin tone color of your skin.

Product code	Size	20 Pieces/Box
CUREDIS-SHER20	19X72 mm	5
	25X72 mm	5
	25Ø Spots	5
	9X38 mm	5




STERILE Or NON STERILE

## Cure - Aid® Clear

**Almost invisible & breathable ideal for working in wet or moist conditions to protect the wound from water.**

Cure - Ad® Clear Plasters with low-adherent absorbent pad placed on transparent plastic washproof material microperforated.

Product code	Size	20 Pieces/Box
CUREDIS-CLER20	19X72 mm	5
	25X72 mm	5
	25Ø Spots	5
	9X38 mm	5



STERILE Or NON STERILE

## Cure - Aid® Sensitive

Extra skin friendly for optimal healing  
suitable for patients with sensitive skin

Cure - Aid® Sensitive Plasters with low-adherent absorbent pad placed on polyester non-woven soft fabric.

Product code	Size	20 Pieces/Box
CUREDIS-SEN20	19X72 mm	5
	25X72 mm	5
	25Ø Spots	5
	9X38 mm	5




## Cure - Aid® Elastic

Ideal for wounds on joints and movable body parts

Cure - Aid® Elastic Plasters with low-adherent absorbent pad placed on fabric and flexible elastic adhesive plaster

Product code	Size	20 Pieces/Box
CUREDIS-ELAS20	19X63 mm	15
	30X64 mm	5



## Cure - Aid® BITE & STING RELIEF patches

Hydrogel bandages due to its high-water content, it provides relief by soothing and cooling and aid patient comfort .

- Immediate cooling relief
- Fast acting cooling effect
- Instant comfort and protection



Product code	Size	10 Pieces/Box
CUREGELADH2045	20X45 mm	10

### Cure - Aid® Cold Sore Patches

The plasters are invisible, makeup or lipstick can be applied over the patch that works on cold sore & acne.

Cure - Aid® Cold Sore Patches with hydrocolloid technology to promote faster healing and leaving no scars on skin.

Product code	Size	24 Pieces/Box
COLLNSPOT24	12mm	9
	10mm	5
	7mm	10

Watch Product Video



### Cure - Aid® Mosquito Repellent

Patch is completely natural plant essence (no harmful chemicals) that repels mosquitoes effectively; it is environmentally friendly and can be used without any side effects.

Attach the adhesive side to your clothes or allocation close to your body such as stockings, shoulder areas, or at edges of the table or bed

**Not on the skin**

Product code	Size	10 Pieces/Box
CUREREP4368	68X43mm	10





**Cure - Aid® Scar Reduction Patch**

Reduces the visibility of raised and coloured scars It is a safe and effective treatment that helps flatten and soften scar tissue.

Cure - Aid® Scar Reduction Patches are a transparent, self-adhesive patch made of polyurethane film with silicone gel sheet made from a clear, soft, medical grade silicone material that sticks to the scar.



Product code	Size	10 Pieces/Box
CURES13872	38X72mm	6


**Cure - Aid® Nose plaster**

Breathe right original nose strips to reduce snoring and relieve nose congestion assists in keep airway open

when applied over the nose and also use with athletes wear nose plaster will increase their air intake and improve their performance.




Product code	Size	10 Pieces/Box
CUREN1967	19X72mm	12



## Cure - Aid® Elastic Finger Plasters

Elastic Finger Plasters It is suitable for covering all types of smaller wounds on fingers

Product code	Size	20 Pieces/Box
CUREDIS-FIN16	H-SHAPE	4
	Fingertip Big	4
	Fingertip small	4
	Finger 2X12	4




STERILE Or NON STERILE

## Cure - Aid® Silver

### Help to prevent infections

The antibacterial silver wound pad is effective against a broad spectrum of bacteria and reduces the risk of infection. Low-adherent absorbent pad placed on skin color plastic micro perforated washproof material.

Product code	Size	12 Pieces/Box
CUREDIS-AG12	19X 72 mm	5
	25 X 72 mm	5
	38 X38 mm	2



STERILE Or NON STERILE

## Cure - Aid® Lipo

**Absorb any exudates to provide moist environment in favour of the healing process.**

Especial gelling non-adherent LIPO-colloid wound pad designed to not stick to any types of wounds placed on polyester non-woven soft fabric

Product code	Size	12 Pieces/Box
CUREDIS-LIPO12	19X 72 mm	5
	25 X 72 mm	5
	38 X38 mm	2



## Cure - Aid® Gentle

**Gentle adhesive allows easy removal with minimal trauma.**

Low-adherent absorbent pad placed on polyester non-woven soft fabric coated with silicone adhesive designed especially for sensitive skin

Product code	Size	12 Pieces/Box
CUREDIS-G12	5	19X 72 mm
	5	25 X 72 mm
	2	38 X38 mm



**Fibrosol® OTC Bandages**

Highly non-adherent absorbent pad made of 100% CMC fibers upon absorption turns to gel creating moist environment without sticking to wound placed on PU waterproof, viral proof & bacterial proof material that protect the wound from external contamination.

Product code	Size	12 Pieces/Box
CUREDIS-FIB12	5	19X 72 mm
	5	25 X 72 mm
	2	38 X38 mm

**Cure - Aid® Algi**

**For bleeding minor cuts, scratches and abrasions**

Stop bleeding haemostatic absorbent primary wound pad placed on polyester non-woven soft fabric

Product code	Size	12 Pieces/Box
CUREDIS-ALGI12	19X 72 mm	5
	25 X 72 mm	5
	38 X38 mm	2

**Cure - Aid® Bamboo** NEW

**The adhesive plaster consists of 100% Biodegradable bamboo**

fabric that is coated with hypoallergenic adhesive and has an absorbent low-adherent pad made of the same material of the Bamboo plaster.

Product code	Size	12 Pieces/Box
CUREDIS-BAM12	19X 72 mm	5
	25 X 72 mm	5
	38 X38 mm	2

Fast Healing



STERILE Or NON STERILE



STERILE Or NON STERILE



STERILE Or NON STERILE

**Cure - Aid® Organic Cotton** NEW

**The adhesive plaster consists of 100% Biodegradable cotton fabric**

that is coated with hypoallergenic adhesive and has an absorbent low-adherent pad made of the same material of the Cotton plaster

Product code	Size	12 Pieces/Box
CUREDIS-COT12	19X 72 mm	5
	25 X 72 mm	5
	38 X38 mm	2

**Cure - Aid® Dressing Strip**

**Classic Rigid Fabric and flexible elastic adhesive plaster durable heavy duty.**

Low-adherent absorbent pad placed on 100% cotton woven Fabric

Product code	Size	10 Pieces/Box
DERMIDIS-CLA10	6cm X 10cm	10



## Cure - Aid® Dressing Strip Washproof

Ideal for working in wet or moist conditions to protect the wound from water.

Low-adherent absorbent pad placed on skin color plastic micro perforated washproof material

Product code	Size	Pieces/Box
DERMIDIS-WP10	6cm X 10cm	10



## Cure - Aid® Dressing Strip Sensitive

Suitable for patients with sensitive skin

Low-adherent absorbent pad placed on Polyester non-woven soft fabric

Product code	Size	Pieces/Box
DERMIDIS-SEN10	6cm X 10cm	10



## Cure - Aid® Dressing Strip Elastic

Ideal for wounds on joints and movable body parts

Low-adherent absorbent pad placed on Flexible Elastic adhesive plaster.

Product code	Size	Pieces/Box
DERMIDIS-ELA10	6cm X 10cm	10



# Detectable Plasters



## Cure - Aid® Blue Metal Detectable Plasters

**Complies with food Hygiene Regulations, waterproof and hypoallergenic**

Supplied with Aluminum Sheet and thin metal strip under the pad which is detectable under the metal detectors placed on PE bright blue color for easier visual detectable

Product code	Size	40 Pieces/Box
CUREDIS-METAL	19x72 mm	10
	19X38 mm	10
	30X 72 mm	10
	38 X72 mm	10



## Cure - Aid® X-ray Detectable Plasters

**Complies with food Hygiene Regulations, waterproof and hypoallergenic**

Supplied with Aluminum Sheet and thin metal strip under the pad which is detectable under the X-ray equipment placed on PE bright blue color for easier visual detectable

Product code	Size	40 Pieces/Box
CUREDIS-XRAY	19x72 mm	10
	19X38 mm	10
	30X 72 mm	10
	38 X72 mm	10



STERILE  OR  NON-STERILE

# Value Packs

## Cure - Aid® Eye Pad

Designed to keep the injured eye safely covered during the healing process  
Island Dressing made of absorbent pad placed centrally on non-woven backing material coated with pressure sensitive adhesive and with opaque black polyurethane layer to protect the eye from light.



STERILE Or NON STERILE

Product	Product code	Size	Pieces/Box
Child	EYEPADN5062	50x62 mm	20
Adult	EYEPADN5882	58x82 mm	20
X-Large (Ocular)	EYEPADN6595	65x95 mm	20

## Cure - Aid® SPOTS

Individually packed first Aid adhesive bandages with absorbent pad placed on skin color plastic micro perforated washproof material.

Box of 100 spot plasters size 25 mm

Product Code	100 Pcs/box
VP-SPOT25	SPOT 25 mm



STERILE Or NON STERILE

## Cure - Aid® eXsorb

Superabsorbent pad placed on skin color plastic micro perforated washproof material.

High absorption capacity absorbs wound exudates and locks it inside the formed gel

Box of 100 plaster assorted that can be customized as your need

Product Code	100 Pcs/box
CUREVALUE-EX100	25x72 mm - 38 X38 mm - 30 X72 mm
	- 38 X72 mm - SPOT 25 mm



STERILE Or NON STERILE



## Cure - Aid® Classic 19x72 mm

Individually packed first Aid adhesive bandages with absorbent pad placed on skin color plastic micro perforated washproof material.

Box of 100 plasters size 19 X 72 mm

Product Code	100 Pcs/box
VP-PE1972	19x72 mm



## Cure - Aid® Value Pack

Individually packed first aid adhesive bandages with absorbent pad placed on different backing materials and different sizes. Box of 100 plaster assorted that can be customized as your need



PACKING MATERIAL Size	PE SKIN COLOR	PE CLEAR	ELASTIC FABRIC	PU	NON – WOVEN	PE BLUE METAL DETECTABLE	BLUE X-RAY	KIDS
19X 38 MM	VP-PE1938	VP-PECL1938	VP-EL1938	VP-PU1938	VP-NW1938	VP-METBLUE1938	VP-XRAY1938	VP-KID1938
19 X 63 MM	VP-PE1963	VP-PECL1963	VP-EL1963	VP-PU1963	VP-NW1963	VP-METBLUE1963	VP-XRAY1963	VP-KID1963
19 X 72 MM	VP-PE1972	VP-PECL1972	VP-EL1972	VP-PU1972	VP-NW1972	VP-METBLUE1972	VP-XRAY1972	VP-KID1972
25 X 72 MM	VP-PE2572	VP-PECL2572	VP-EL2572	VP-PU2572	VP-NW2572	VP-METBLUE2572	VP-XRAY2572	VP-KID2572
38 X 38 MM	VP-PE3838	VP-PECL3838	VP-EL3838	VP-PU3838	VP-NW3838	VP-METBLUE3838	VP-XRAY3838	VP-KID3838
ANCHOR	VP-PE3072	VP-PECL3072	VP-EL3072	VP-PU3072	VP-NW3072	VP-METBLUE3072	VP-XRAY3072	VP-KID3072
BUTTERFLY	VP-PE4563	VP-PECL4563	VP-EL4563	VP-PU4563	VP-NW4563	VP-METBLUE4563	VP-XRAY4563	VP-KID4563
50 X 72 MM	VP-PE5072	VP-PECL5072	VP-EL5072	VP-PU5072	VP-NW5072	VP-METBLUE5072	VP-XRAY5072	VP-KID5072
38 X 72 MM	VP-PE3872	VP-PECL3872	VP-EL3872	VP-PU3872	VP-NW3872	VP-METBLUE3872	VP-XRAY3872	VP-KID3872