

Cure - Aid[®]



Description:

Cure-Aid[®] plasters and dressing strip: adhesive plasters that consists of an absorbent low-adherent pad placed on different backing materials coated with pressure sensitive adhesive.

Cure-Aid[®] Blue Metal & X-ray Detectable plasters and dressing strip: supplied with an Aluminum sheet and thin metal strip under the pad which is detectable in most metal detectors and X-ray equipment. Bright blue color for easier visual detection.

Cure-Aid[®] Eye pad: Island dressing made of absorbent pad, placed centrally on non-woven backing material coated with pressure sensitive adhesive and with opaque black polyurethane layer to protect the eye from light. Designed to keep the injured eye safely covered during the healing process.


Cure-Aid[®] is available with different backing materials:

- 1- **Polyethylene Foam:** semipermeable, waterproof and with small open cells that can hold fluid.
- 2- **Polyethylene film:** perforated film supplied with multiple colors e.g. transparent, white, skin color.
- 3- **Printed films:** waterproof and dirt-repellent kids plaster with fun motifs.
- 4- **Polyurethane film:** waterproof, bacteria-proof and virus-proof plasters that protect the wound from external contamination.
- 5- **Non-woven:** soft fabric suitable for patients with sensitive skin.
- 6- **Non-woven Ultra:** additional upper polyurethane layer laminated to the non-woven backing, protects the wound from external contamination and suitable for sensitive skin.
- 7- **Universal:** additional upper polyethylene layer laminated to the non-woven backing, protects the wound from external moisture and suitable for sensitive skin.
- 8- **Thermoplastic Polyurethane:** TPU makes the plaster breathable and flexible.
- 9- **100% Cotton Fabric:** Low irritant Rigid Fabric and Flexible Elastic adhesive plaster. Supplied in special shapes designed for fingertips and joints. Also supplied as Side Way Elastic Fabric (transversely stretch-able) that adapts well even to highly contoured parts of the body. Stretch is very well tolerated even by patients with sensitive skin.
- 10- **Paper Non-Woven:** medical paper backing suitable for patients with sensitive skin.
- 11- **PU Tough:** scrim net substrate laminated with Polyurethane film which makes the plaster rigid.

Cure-Aid[®] Plasters and dressing strip: product line is available in sterile and non-sterile variants


Cure-Aid[®] is available with different applications systems:

- 1- Wave-cut application system.




Cure - Aid[®]

2- Frame Style application system.



3-Cure-aid easy open.



Features:

Cure-Aid[®] plasters and dressing strip:

- Breathable with high permeability to water vapor and oxygen.
- Low-adherent absorbent pad that doesn't stick to the wound.
- Highly conformable and flexible plasters suitable for curves and contours of body.
- Easy to apply and with rounded edges to prevent dressing from curling up.

Cure-Aid[®] Blue Metal & X-ray Detectable plasters and dressing strip:

- A necessity in food-processing plants where a foreign object in finished product can be a disaster
- Waterproof.
- Visually Detectable.
- X-Ray and metal detectable.

Cure-Aid[®] Eye pad:

- Protects the eyes from dust and grime.
- Soft and comfortable absorbent pad.
- Strong elastic bandage secures pad firmly in place.

Intended Purpose:

Cure-Aid[®] plasters and dressing strip:

Intended for the protection of minor acute wounds such as cuts, skin scratches, abrasions, lacerations and after blood drawing or vaccination.

Cure-Aid[®] Blue Metal & X-ray Detectable plasters and dressing strip:

Ideal to cover minor wounds, abrasions, superficial wounds and skin cuts in food processing factories, restaurants,...etc.

Cure-Aid[®] Eye pad:

Intended for patients with minor eye injuries, foreign objects embedded in their eyes and with lazy eye disorder. Besides, the pad can be used after eye trauma e.g. corneal abrasion and may provide initial protection from possible risk of infection.

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263

- Kafr-El-Zayat P.O.Box 38 Egypt

Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com

 **Pharmaplast**

www.pharmaplast-online.com



Embossed Pu Blue



Embossed Pu Beige



Jungle




Neon Orange



Neon Green



Neon Rose



Premium rigid Fabric



Skin Colour Sensitive



Foam



Plastic Washproof



Blue Washproof Detectable



Dark brown Washproof



Plastic Clear



Blue Elastic Fabric Detectable




Beige Nonwoven



Pu Film



Elastic Fabric



Nonwoven




Kids Printed



Rigid Fabric



Sensitive




Universal



Premium Elastic Fabric



Children animals plaster



Tough Plaster



CUREN3055 30 x 55 mm



CUREN3072 30 x 72 mm



Eye Pad adult
EYEPADN5882 58 x 82 mm



Eye Pad child
EYEPADN5062 50 x 62 mm








Oval adult
CURENKNEE 64 x 110 mm








Oval child
CUREN5284 52 x 84 mm





























Cure - Aid[®]

Cure-Aid [®] Adhesive plaster	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
CURE938 	9 X 38 mm	0.35 X 1.5 in	100	60	56 X 29 X 29 cm
CURE1657 	16 X 57 mm	0.63 X 2.24 in	100	60	56 X 29 X 29 cm
CURE1938 	19 X 38 mm	0.75 X 1.5 in	100	60	56 X 29 X 29 cm
CURE1963 	19 X 63 mm	0.75 X 2.48 in	100	60	56 X 29 X 29 cm
CURE1964 	19 X 64 mm	0.75 X 2.52 in	100	60	56 X 29 X 29 cm
CURE1972 	19 X 72 mm	0.75 X 2.83 in	100	60	67.5 X 30 X 33 cm
CURE1976 	19 X 76 mm	0.75 X 2.99 in	100	60	56 X 29 X 29 cm
CURE2572 	25 X 72 mm	0.98 X 2.83 in	100	60	56 X 29 X 29 cm
CURE2576 	25 X 76 mm	0.98 X 2.99 in	100	60	56 X 29 X 29 cm
CURE3055 	30 X 55 mm	1.18 X 2.17 in	100	60	56 X 29 X 29 cm
CURE3064 	30 X 64 mm	1.18 X 2.52 in	100	60	56 X 29 X 29 cm
CURE3072 	30 X 72 mm	1.18 x 2.83 in	50	60	56 X 29 X 29 cm
CUREANCHOR 	30 X 72 mm Anchor	1.18 x 2.84 in	50	60	56 X 29 X 29 cm
CURE3838 	38 X 38 mm	1.5 X 1.5 in	100	60	56 X 29 X 29 cm
CURE3872 	38 X 72 mm	1.5 X 2.83 in	50	60	56 X 29 X 29 cm
CURE5072 	50 X 72 mm	1.97 X 2.83 in	50	60	56 X 29 X 29 cm
CURE50100 	50 X 100 mm	2 X 4 in	50	60	56 X 29 X 29 cm
CURE5060 	50 X 60 mm Fingertip	1.97 X 2.36 in	50	60	56 X 29 X 29 cm
CURE4563 	45 X 63 mm Butterfly	1.77 X 2.48 in	50	60	56 X 29 X 29 cm
CURE25 	SPOT 25 mm	0.98 in	100	192	58 X 33 X 30 cm
CURE22 	SPOT 22 mm	0.8 in	100	192	58 X 33 X 30 cm
CURE2525 	25 X 25 mm	0.98 X 0.98 in	100	60	56 X 29 X 29 cm
EYEPAD5062 	50x62 mm Eyepad child	Eyepad child	20	40	47 X 23 X 20 cm
EYEPAD5882 	58x82 mm Eyepad Adult	Eyepad Adult	20	40	47 X 23 X 20 cm
EYEPAD6595 	65x95 mm Ocular	Ocular	20	40	47 X 23 X 20 cm
CUREFEET 	Feet	Feet	10	48	
CURE2235 	22 x 35 mm	0.9 x 1.4 in	100	60	56 X 29 X 29 cm
CURE2665 	26 x 65 mm	1.5 x 2.6 in	100	60	56 X 29 X 29 cm
CURE3070 	30 x 70 mm	1.18 x 2.8 in	100	60	56 X 29 X 29 cm

Cure - Aid®

Cure-Aid® Adhesive plaster	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
CURE64110 	64 X 110 mm OVALADULT	2.5 x 4.3 in	50	60	56 X 29 X 29 cm
CURE5284 	52 X 84 mm OVALCHILD	2.04 x 3.3 in	50	60	56 X 29 X 29 cm
CURE20120 	20 x 120 mm	0.8 x 4.7 in	100	60	56 X 29 X 29 cm
CURE30120 	30 x 120 mm	1.2 x 4.7 in	100	60	56 X 29 X 29 cm
CURE20180 	20 x 180 mm	0.8 x 7 in	100	60	56 X 29 X 29 cm
CURE30180 	30 x 180 mm	1.2 x 7 in	100	60	56 X 29 X 29 cm
CURE2582 	25 x 82 mm	1 x 3.2 in	100	60	56 X 29 X 29 cm
CURE2956 	29 x 56 mm	1.1 x 2.2 in	100	60	56 X 29 X 29 cm
CURECHARACTERS 	19x63 mm	0.75 X 2.48 in	100	60	56 X 29 X 29 cm
CURE1955 	19 x 55 mm	0.7 x 2.2 in	100	60	56 X 29 X 29 cm
CURE4550 	45 x 50 mm Butterfly	1.8 x 2.2 in	50	60	56 X 29 X 29 cm
CURE1966 	19x66 mm	0.75 x 2.6 in	100	60	
CURE1045 	10x45 cm	0.4 x 1.8 in	100	60	
CUREHEART 	Heart		100	60	
CURELIPS 	Lips		100	60	
CUREKNEE 	knee		50	60	
Cure-Aid® Dermi Plaster	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
DERMI60100 	60 x 100 mm	2.4 x 4 in	10	48	
DERMI601000 	60 x 1000 mm	2.4 x 39.4 in	1	48	
DERMI80100 	80 X 100 mm	3.1 x 4 in	10	120	
DERMI801000 	80 x 1000 mm	3.1 x 39.4 in	1	120	
DERMI40100 	40 x 100 mm	1.57 x 4 in	10	120	
DERMI60200 	60 x 200 mm	2.4 x 7.87 in	10	120	
DERMI2560 	25 x 60 mm	0.98 x 2.4 in	100	60	56 X 29 X 29 cm
DERMI401000 	40 x 1000 mm	1.57 x 39.4 in	1	120	
DERMI40500 	40 x 500 mm	1.57 x 19.69 in	1	120	
DERMI60500 	60 x 500 mm	2.4 x 19.69 in	1	120	
DERMI80500 	80 x 500 mm	3.15 x 19.69 in	1	120	
DERMI405000 	40 x 5000 mm	1.57 x 196.85 in	1	40	
DERMI605000 	60 x 5000 mm	2.4 x 196.85 in	1	40	

Cure - Aid®

Cure-Aid® Dermi plaster	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
DERMI805000 	80 x 5000 mm	3.1 x 196.85 in	1	40	
DERMI20120 	20 x 120 mm	0.8 x 4.7 in	100	60	56 X 29 X 29 cm
DERMI30120 	30 x 120 mm	1.18 x 4.7 in	100	60	56 X 29 X 29 cm
DERMI20180 	20 x 180 mm	0.8 x 7 in	100	60	56 X 29 X 29 cm
DERMI30180 	30 x 180 mm	1.18 x 7 in	100	60	56 X 29 X 29 cm
DERMI20120IW 	20 x 120 mm	0.8 x 4.7 in	100	60	56 X 29 X 29 cm
DERMI30120IW 	30 x 120 mm	1.18 x 4.7 in	100	60	56 X 29 X 29 cm
DERMI20180IW 	20 x 180 mm	0.8 x 7 in	100	60	56 X 29 X 29 cm
DERMI30180IW 	30 x 180 mm	1.18 x 7 in	100	60	56 X 29 X 29 cm
DERMI2040 	20x40 mm	0.8 x 7 in	100	60	56 X 29 X 29 cm
DERMI2060 	20x60 mm	0.8 x 2.4 in	100	60	56 X 29 X 29 cm
DERMI70100 	70x100 mm	2.8 x 4 in	50	60	
Cure-Aid® Plaster Non-Sterile	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
CUREN938 	9 X 38 mm	0.35 X 1.5 in	100	60	56 X 29 X 29 cm
CUREN1657 	16 X 57 mm	0.63 X 2.24 in	100	60	56 X 29 X 29 cm
CUREN1938 	19 X 38 mm	0.75 X 1.5 in	100	60	56 X 29 X 29 cm
CUREN1963 	19 X 63 mm	0.75 X 2.48 in	100	60	56 X 29 X 29 cm
CUREN1964 	19 X 64 mm	0.75 X 2.52 in	100	60	56 X 29 X 29 cm
CUREN1972 	19 X 72 mm	0.75 X 2.83 in	100	60	67.5 X 30 X 33 cm
CUREN1976 	19 X 76 mm	0.75 X 2.99 in	100	60	56 X 29 X 29 cm
CUREN2572 	25 X 72 mm	0.98 X 2.83 in	100	60	56 X 29 X 29 cm
CUREN2576 	25 X 76 mm	0.98 X 2.99 in	100	60	56 X 29 X 29 cm
CUREN3055 	30 X 55 mm	1.18 X 2.17 in	100	60	56 X 29 X 29 cm
CUREN3064 	30 X 64 mm	1.18 X 2.52 in	100	60	56 X 29 X 29 cm
CUREN3072 	30 X 72 mm	1.18 x 2.83 in	50	60	56 X 29 X 29 cm
CURENANCHOR 	30 X 72 mm Anchor	1.18 x 2.84 in	50	60	56 X 29 X 29 cm
CUREN3838 	38 X 38 mm	1.5 X 1.5 in	100	60	56 X 29 X 29 cm
CUREN3872 	38 X 72 mm	1.5 X 2.83 in	50	60	56 X 29 X 29 cm
CUREN5072 	50 X 72 mm	1.97 X 2.83 in	50	60	56 X 29 X 29 cm

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263

- Kafr-El-Zayat P.O.Box 38 Egypt





























Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com

 **Pharmaplast**

www.pharmaplast-online.com

Cure - Aid®

Cure-Aid® Plaster Non-Sterile	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
CUREN50100 	50 X 100 mm	2 X 4 in	50	60	56 X 29 X 29 cm
CUREN5060 	50 X 60 mm Fingertip	1.97 X 2.36 in	50	60	56 X 29 X 29 cm
CUREN4563 	45 X 63 mm Butterfly	1.77 X 2.48 in	50	60	56 X 29 X 29 cm
CUREN25 	SPOT 25 mm	0.98 in	100	192	58 X 33 X 30 cm
CUREN22 	SPOT 22 mm	0.8 in	100	192	58 X 33 X 30 cm
CUREN2525 	25 X 25 mm	0.98 X 0.98 in	100	60	56 X 29 X 29 cm
EYEPADN5062 	50x62 mm Eyepad child	Eyepad child	20	40	47 X 23 X 20 cm
EYEPADN5882 	58x82 mm Eyepad Adult	Eyepad Adult	20	40	47 X 23 X 20 cm
EYEPADN6595 	65x95 mm Ocular	Ocular	20	40	47 X 23 X 20 cm
CURENFEET 	Feet	Feet	10	48	
CUREN2235 	22 x 35 mm	0.9 x 1.4 in	100	60	56 X 29 X 29 cm
CUREN2665 	26 x 65 mm	1.5 x 2.6 in	100	60	56 X 29 X 29 cm
CUREN3070 	30 x 70 mm	1.18 x 2.8 in	100	60	56 X 29 X 29 cm
CUREN64110 	64 X 110 mm OVALADULT	2.5 x 4.3 in	50	60	56 X 29 X 29 cm
CUREN5284 	52 X 84 mm OVALCHILD	2.04 x 3.3 in	50	60	56 X 29 X 29 cm
CUREN20120 	20 x 120 mm	0.8 x 4.7 in	100	60	56 X 29 X 29 cm
CUREN30120 	30 x 120 mm	1.2 x 4.7 in	100	60	56 X 29 X 29 cm
CUREN20180 	20 x 180 mm	0.8 x 7 in	100	60	56 X 29 X 29 cm
CUREN30180 	30 x 180 mm	1.2 x 7 in	100	60	56 X 29 X 29 cm
CUREN2582 	25 x 82 mm	1 x 3.2 in	100	60	56 X 29 X 29 cm
CUREN2956 	29 x 56 mm	1.1 x 2.2 in	100	60	56 X 29 X 29 cm
CURENCHARACTERS 	CHARACTERS	0.75 X 2.48 in	100	60	56 X 29 X 29 cm
CUREN1955 	19 x 55 mm	0.7 x 2.2 in	100	60	56 X 29 X 29 cm
CUREN4550 	45 x 50 mm Butterfly	1.8 x 2.2 in	50	60	56 X 29 X 29 cm
CUREN3868 	38 x68 mm Butterfly	1.5 x 2.7 in	50	60	56 X 29 X 29 cm
CUREN60100 	60 x 100 mm	2.4 x 4 in	10	48	56 X 29 X 29 cm
CUREN601000 	60 x 1000 mm	2.4 x 39.4 in	1	48	56 X 29 X 29 cm
CUREN3272 	32x72 mm	1.26 x 2.8 in	50	60	56 X 29 X 29 cm
CUREN1966 	19x66 mm	0.75 x 2.59 in	100	60	56 X 29 X 29 cm

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263

- Kafr-El-Zayat P.O.Box 38 Egypt





























Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com

 **Pharmaplast**

www.pharmaplast-online.com

Cure - Aid®

Cure-Aid® Plaster Non-Sterile	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
CUREN1045 	10x45 mm	0.39 x 1.77 in	100	60	56 X 29 X 29 cm
CURENHEART 	Heart		100	60	56 X 29 X 29 cm
CURENLIPS 	Lips		100	60	56 X 29 X 29 cm
CURENKNEE 	64x110 mm	2.5 x 4.33 in	50	60	
Cure-Aid® Dermi Plaster Non-Sterile	size	Size (Inch)	pieces / box	Boxes/Case	Case Dimensions
DERMIN2040 	20x40 mm	0.8 x 1.6 in	100	60	56 X 29 X 29 cm
DERMIN2060 	20x60 mm	0.8x 2.4 in	100	60	56 X 29 X 29 cm
DERMIN20120 	20x120 mm	0.8 x 4.7 in	100	60	56 X 29 X 29 cm
DERMIN30120 	30 x 120 mm	1.2 x 4.7 in	100	60	56 X 29 X 29 cm
DERMIN20180 	20 x 180 mm	0.8 x 7 in	100	60	56 X 29 X 29 cm
DERMIN30180 	30 x 180 mm	1.2 x 7 in	100	60	56 X 29 X 29 cm
DERMIN40100 	40x100 mm	1.6 x 4 in	10	120	
DERMIN60100 	60x100 mm	2.4x4 in	10	120	
DERMIN601000 	60x1000 mm	2.4x39.4 in	1	120	
DERMIN70100 	70x100 mm	2.7 x 4 in	10	120	
DERMIN80100 	80 X 100 mm	3.1 x 4 in	10	120	
DERMIN801000 	80 X 1000 mm	3.1 x 39.4 in	1	120	
DERMIN605000 	60X5000 mm	2.4 x 197 in	1	40	
DERMIN805000 	80X5000 mm	3.1 x 197 in	1	120	
DERMIN60500 	60x500 mm	2.4 x 19.69 in	1	120	
DERMIN80500 	80x500 mm	3 x 19.69 in	1	40	
DERMIN60200 	60x200 mm	2.4 x 8.7 in	10	120	
DERMIN2560 	25x60 mm	0.98 x 2.4 in	100	60	
DERMIN40500 	40x500 mm	1.8 x 19.69 in	1	120	
DERMIN405000 	40x5000 mm	1.8 x 196.9 in	1	40	
DERMIN80750 	80x750 mm	3.14 x 29.52 in	1	120	
DERMIN603000 	60x3000 mm	2.36 x 118.11 in	1	48	
DERMIN604000 	60x4000 mm	2.36 x 157.48 in	1	48	
DERMIN20120IW 	20 x 120 mm individually wrapped	0.8 x 4.7 in	100	60	
DERMIN30120IW	30 x 120 mm individually wrapped	1.18 x 7 in	100	60	
DERMIN20180IW	20 x 180 mm individually wrapped	0.8 x 7 in	100	60	
DERMIN30180IW	30 x 180 mm individually wrapped	1.18 x 7 in	100	60	

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263


- Kafr-El-Zayat P.O.Box 38 Egypt

Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com

 **Pharmaplast**

www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Cure-Aid						
Issue date : 19-10-2012		Prepared by : Marina Medhat						
Colors : 		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	05-05-2015	new text by Jostina Ramsis	Essential requirement					
3	13-11-2021	updating codes	to be aligned with new coding system					
4	27-01-2023	Adding new Codes	to be aligned with new coding system					