

Difa IV[®] CHG Catheter Securement



www.zaditibbigerecler.com

Difa IV[®] CHG Catheter Securement



Description:

Difa IV[®] Catheter securement dressings are a vast range of products that use the **Difa IV[®]** gel technology with chlorhexidine to secure catheters, hoses and pipes. **Difa IV[®]** gel is a water clear hydrogel, made of ionic polymer, chlorhexidine and glycerol it is placed as an island pad on a PU film backing coated with hypoallergenic adhesive.

Difa IV[®] provides soft contact layer with high absorption capacity.

Features:

- Highly absorbent gel layer.
- Invitro Studies show bacteriostatic and fungistatic effect up to 7 days.
- Tackiness: adheres to skin, so it can be applied without fixation tapes of secondary film dressings.
- Easy removal: it can be removed without pain, trauma or residues on skin.

Benefits:

1. **Difa IV[®] Catheter securement** dressings are comprehensively designed to deliver exceptional patient care through advanced catheter securement, superior wear time, gentle removal, and excellent absorption.
2. **Difa IV[®] Catheter securement** dressings help secure and preserve the integrity of devices, minimize movement, prevent catheter dislodgement.
3. Invitro tests have shown that **Difa IV[®]** exhibits bacteriostatic and fungistatic effect.
4. **Difa IV[®] Catheter securement** dressings are ideal for use during catheter application for prevention of intra-vascular catheter-associated infections

Indications:

• Difa IV dressings:

1. Cannula and Catheters fixation.

Contraindications:

None










Difa IV[®] CHG Catheter Securement

How Supplied:

Difa IV[®] Catheter securement dressings are supplied as sterile, individually wrapped in medical sealed paper.

Difa IV[®] Catheter securement dressings are sterilized by EO gas.

Sizes

Difa PU IV CHG	Size	Pieces/ box	Boxes / Case
Difa IV CHG 76	7 cm x 6 cm	100	16
Difa IV CHG 77	7 cm x 7 cm	100	16
Difa IV CHG 87	8 cm x 7 cm	100	16
Difa IV CHG 810	8 cm x 10 cm	100	16
Difa IV CHG 2810	8.25 cm x 10 cm	100	16
Difa IV CHG 856	8.5 cm x 6 cm	100	16
Difa IV CHG 857	8.5 cm x 7 cm	100	16
Difa IV CHG 8510	8.5 cm x 10 cm	100	16
Difa IV CHG 97	9 cm x 7 cm	100	16
Difa IV CHG 910	9 cm x 10 cm	100	16
Difa IV CHG 1012	10 cm x 12 cm	100	16
Difa IV CHG 1014	10 m x 14 cm	100	16
Difa IV CHG 68	6 cm x 8 cm	100	16
Difa PU IV FS CHG	Size	Pieces/ box	Boxes / Case
Difa IV CHG 79 	9 cm x 7 cm	100	16
Difa IV CHG 85115 	8.5 cm x 11.5 cm	100	16
Difa IV CHG 5057 	5 cm x 5.7 cm	100	16
Difa IV CHG 67 	6 cm x 7 cm	100	16
Difa IV CHG 657 ADV 	6.5 cm x 7 cm	25	24
Difa IV CHG 79 ADV 	7 cm x 9 cm	25	24
Difa IVCHG85115 ADV 	8.5 cm x 11.5 cm	25	24
Difa IVCHG1012 ADV 	10 cm x 12 cm	25	24
Difa IVCHG10155 ADV 	10 cm x 15.5 cm	25	24

Zadi Tıbbi Gereçler San. ve Tic. A.Ş
İOSB Mah. İmsan Sanayi Sit. B Blok No:13
Küçükçekmece / İSTANBUL
Tel: +90 212 495 73 03 Fax: +90 212 495 73 04
Email: import@zaditibbigerecler.com

www.zaditibbigerecler.com