

Gentle Post-operative Dressing & OTC bandages







Post-operative silicone dressings:

Pharmapore® Gentle: Post-Operative dressing made of soft and conformable nonwoven fixative layer, coated with silicone adhesive and absorbent wound pad laid at the center of the nonwoven layer. The wound contact side is protected with two release liner flaps. The dressing is individually packed in heat seal pouches and sterile by EO gas.

Pharmapore® PU Gentle: Post-operative dressing made of waterproof, virus proof and bacteria proof polyurethane film that allows showering while the dressing is in place. The polyurethane film comes with top backing paper to allow application of the film without wrinkling.

Pharmapore® PU Frame Style Gentle: Same as **Pharmapore PU Gentle** but with frame delivery system to allow easy application.

Pharmapore® Ultra Gentle: Post-operative Dressing made of polyurethane film layer laminated on top of the nonwoven layer, coated with Silicone adhesive and absorbent wound pad placed at the center. It combines the softness of Pharmapore Gentle and the waterproofness of **Pharmapore® PU Gentle**.



Gentle Post-operative Dressing & OTC bandages

Silicone Film Dressings:

Protectfilm® Gentle: Transparent, waterproof, virus proof and bacteria proof film coated with silicone adhesive.






Protectfilm® Gentle Frame style: Transparent, waterproof, virus proof and bacteria proof film coated with silicone adhesive, with frame delivery system for easy application.

IV silicone dressings:

Pharmapore® IV Gentle: Soft conformable nonwoven layer, coated with silicone adhesive and absorbent wound pad placed at the center, having U shaped slit. The structure allows easy fixation of catheters and cannula.

Pharmapore® PU IV Gentle: Waterproof, Virus proof and bacteria proof polyurethane film coated with silicone adhesive and absorbent wound pad placed at the center, having U shaped slit. The structure allows easy fixation of catheters and cannula.


Pharmapore® PU IV Gentle FS: Same as **Pharmapore® PU IV Gentle**, but with frame delivery system for easy application.



Gentle Post-operative Dressing & OTC bandages

Silicone adhesive OTC bandages:

Cure-Aid® Gentle: Adhesive plasters for the retail market made of different types of substrates e.g. polyurethane film, polyethylene film, elastic fabric, nonwoven,...etc., coated with silicone adhesive and absorbent wound pad at the center of the substrate. The Cure-Aid Gentle is individually packed and sterilized by EO gas.



Indication:

Post-operative silicone dressings:

Post-operative silicone dressings are designed for exuding wounds. They are intended for acute wounds, such as surgical wounds, cuts and abrasions. It is optimized for post-op use and blood absorption.

IV silicone dressings

IV silicone dressings are designed for cannula and catheter fixation, allowing atraumatic removal or dressing change without exposing the patient to additional pain.

Silicone Film Dressings

Silicone Film Dressings are designed for a wide range of superficial wounds such as pressure ulcers category 1 and 2, superficial skin injuries and superficial burns. Silicone Film Dressings protects fragile and sensitive skin.

Silicone adhesive plasters:

Indications: minor cuts, scratches and abrasions

Contraindications:

· Dressings should not be used on people who are allergic to the dressing or any of its components.

Precautions:

- Dressings are not intended for use within internal body cavities.
- If infection develops during the use of the dressing, clinically indicated antibiotic therapy should be initiated.
- Control of blood glucose should be monitored for diabetic foot ulcers.
- Delicate newly formed blood vessels could produce blood stained wound fluid.
- Sterility is guaranteed unless pouch is damaged.

Instructions for use:

- Initial disinfection of hands before Gentle dressing application.
- Cleanse the wound.
- The dressing should overlap at least 1 cm onto the skin surrounding the wound.

Removal:

- Dressing may remain on the wound for up to 7 days.
- Dressing may be removed when wound is about to close.
- The dressing should be removed when clinically indicated.
- Irrigate with saline if residual dressing is left on the wound.

Features:

- Soft breathable cover
- Highly absorbent and conformable pad

Gentle Post-operative Dressing & OTC bandages

- Low allergy adhesive
- Safe and secure fixation
- Gentle adhesive allows easy removal with minimal trauma

Sizes:

Pharmapore PU IV Gentle Frame Style (with Pad)	Pharmapore PU IV Gentle Frame Style	Pharmapore PU IV Gentle	Pharmapore IV Gentle	size	Pieces / box	Boxes / Case
IVFSGP6070	IVFSG6070	IVPUG6070	IVNWG6070	6X7CM	100	16
IVFSGP7070	IVFSG7070	IVPUG7070	IVNWG7070	7X7CM	100	16
IVFSGP7080	IVFSG7080	IVPUG7080	IVNWG7080	7X8CM	100	16
IVFSGP80100	IVFSG80100	IVPUG80100	IVNWG80100	8X10CM	100	16
IVFSGP82100	IVFSG82100	IVPUG82100	IVNWG82100	8.25X10CM	100	16
IVFSGP6085	IVFSG6085	IVPUG6085	IVNWG6085	6X8.5CM	100	16
IVFSGP7085	IVFSG7085	IVPUG7085	IVNWG7085	7X8.5CM	100	16
IVFSGP85100	IVFSG85100	IVPUG85100	IVNWG85100	8.5X10CM	100	16
IVFSGP85115	IVFSG85115	IVPUG85115	IVNWG85115	8.5X11.5CM	100	16
IVFSGP7090	IVFSG7090	IVPUG7090	IVNWG7090	7X9CM	100	16
IVFSGP90100	IVFSG90100	IVPUG90100	IVNWG90100	9X10CM	100	16
IVFSGP100120	IVFSG100120	IVPUG100120	IVNWG100120	10X12CM	50	16
IVFSGP100140	IVFSG100140	IVPUG100140	IVNWG100140	10X14CM	50	16
IVFSGP6080	IVFSG6080	IVPUG6080	IVNWG6080	6X8CM	100	16
IVFSGP5057	IVFSG5057	IVPUG5057	IVNWG5057	5X5.7CM	100	16
	IVFSG100120ADV		IVNWG100155	10x12CM	50	24
	IVFSG100155ADV		IVNWG5055	10x15.5CM	50	24
	IVFSG6570ADV		IVNWG6570	6.5x7CM	50	24
	IVFSG7085ADV			7x8.5CM	50	24
	IVFSG85115ADV			8.5x11.5CM	50	24
Pharmapore Gentle	Pharmapore Ultra Gentle	Pharmapore PU Gentle	Pharmapore PU Gentle FS	Size	Pieces / box	Boxes / Case
POREG5070	POREG5070	POREPUG5070	POREPUGFS5070	5cmx7cm	10	48
POREG6070	POREG6070	POREPUG6070	POREPUGFS6070	6cmx7cm	10	48
POREG6080	POREG6080	POREPUG6080	POREPUGFS6080	6cmx8cm	200	8
					10	48
POREG5090	POREG5090	POREPUG5090	POREPUGFS5090	5cmx9cm	10	48
POREG6083	POREG6083	POREPUG6083	POREPUGFS6083	6cmx8.3cm	10	48
POREG50100	POREG50100	POREPUG50100	POREPUGFS50100	5cmx10cm	10	48
POREG6082	POREG6082	POREPUG6082	POREPUGFS6082	6cmx8.25cm	10	48
POREG6085	POREG6085	POREPUG6085	POREPUGFS6085	6cmx8.5cm	10	48
POREG7080	POREG7080	POREPUG7080	POREPUGFS7080	7cmx8cm	10	48
POREG6580	POREG6580	POREPUG6580	POREPUGFS6580	6.5cmx8cm	10	48
POREG7575	POREG7575	POREPUG7575	POREPUGFS7575	7.5cmx7.5cm	10	48
POREG6090	POREG6090	POREPUG6090	POREPUGFS6090	6cmx9cm	10	48
POREG60100	POREG60100	POREPUG60100	POREPUGFS60100	6cmx10cm	10	48
POREG5082	POREG5082	POREPUG5082	POREPUGFS5082	5cmx8.25cm	10	48
POREG75100	POREG75100	POREPUG75100	POREPUGFS75100	7.5cmx10cm	10	48
POREG80100	POREG80100	POREPUG80100	POREPUGFS80100	8cmx10cm	10	48
POREG82100	POREG82100	POREPUG82100	POREPUGFS82100	8.25cmx10cm	10	48
POREG100100	POREG100100	POREPUG100100	POREPUGFS100100	10cmx10cm	10	48
POREG70120	POREG70120	POREPUG70120	POREPUGFS70120	7cmx12cm	10	48
POREG80120	POREG80120	POREPUG80120	POREPUGFS80120	8cmx12cm	10	48
POREG100120	POREG100120	POREPUG100120	POREPUGFS100120	10cmx12cm	10	48
POREG83120	POREG83120	POREPUG83120	POREPUGFS83120	8.3cmx12cm	10	48
POREG90150	POREG90150	POREPUG90150	POREPUGFS90150	9cmx15cm	10	48
POREG80150	POREG80150	POREPUG80150	POREPUGFS80150	8cmx15cm	10	48
POREG82150	POREG82150	POREPUG82150	POREPUGFS82150	8.25cmx15cm	10	48
POREG100140	POREG100140	POREPUG100140	POREPUGFS100140	10cmx14cm	10	48
POREG100150	POREG100150	POREPUG100150	POREPUGFS100150	10cmx15cm	10	48
POREG83180	POREG83180	POREPUG83180	POREPUGFS83180	8.3cmx18cm	10	48
POREG82200	POREG82200	POREPUG82200	POREPUGFS82200	8.25cmx20cm	10	24
POREG90200	POREG90200	POREPUG90200	POREPUGFS90200	9cmx20cm	10	24
POREG100200	POREG100200	POREPUG100200	POREPUGFS100200	10cmx20cm	10	24
POREG83240	POREG83240	POREPUG83240	POREPUGFS83240	8.3cmx24cm	10	24
POREG100120	POREG100120	POREPUG100120	POREPUGFS100120	10cmx12cm	10	48

Gentle Post-operative Dressing & OTC bandages

Sizes:

Pharmapore Gentle	Pharmapore Ultra Gentle	Pharmapore PU Gentle	Pharmapore PU Gentle FS	size	Pieces / box	Boxes / Case
POREG82250	POREG82250	POREPUG82250	POREPUGFS82250	8.25cmx25cm	10	24
POREG90250	POREG90250	POREPUG90250	POREPUGFS90250	9cmx25cm	10	24
POREG100250	POREG100250	POREPUG100250	POREPUGFS100250	10cmx25cm	10	24
POREG82300	POREG82300	POREPUG82300	POREPUGFS82300	8.25cmx30cm	10	24
POREG90300	POREG90300	POREPUG90300	POREPUGFS90300	9cmx30cm	10	24
POREG100300	POREG100300	POREPUG100300	POREPUGFS100300	10cmx30cm	10	24
POREG90350	POREG90350	POREPUG90350	POREPUGFS90350	9cmx35cm	10	24
POREG82350	POREG82350	POREPUG82350	POREPUGFS82350	8.25cmx35cm	10	24
POREG100350	POREG100350	POREPUG100350	POREPUGFS100350	10cmx35cm	10	24
POREG83415	POREG83415	POREPUG83415	POREPUGFS83415	8.3cmx41.5cm	10	24
POREG100160	POREG100160	POREPUG100160	POREPUGFS100160	10cmx16cm	10	48
POREG100210	POREG100210	POREPUG100210	POREPUGFS100210	10cmx21cm	10	24
POREG100260	POREG100260	POREPUG100260	POREPUGFS100260	10cmx26cm	10	24
POREG100400	POREG100400	POREPUG100400	POREPUGFS100400	10cmx40cm	10	24
POREG120160	POREG120160	POREPUG120160	POREPUGFS120160	12cmx16cm	10	48
POREG150150	POREG150150	POREPUG150150	POREPUGFS150150	15cmx15cm	10	24
POREG200200	POREG200200	POREPUG200200	POREPUGFS200200	20cmx20cm	10	24
POREG2572	POREG2572	POREPUG2572	POREPUGFS2572	2.5cmx7.2cm	400	8
POREG38100	POREG38100	POREPUG38100	POREPUGFS38100	8.3cmx10cm	10	48
POREG5072	POREG5072	POREPUG5072	POREPUGFS5072	5cmx7.2cm	10	48
POREG5075	POREG5075	POREPUG5075	POREPUGFS5075	5cmx7.5cm	10	48
POREG70100	POREG70100	POREPUG70100	POREPUGFS70100	7cmx10cm	10	48
POREG6075	POREG6075	POREPUG6075	POREPUGFS6075	6cmx7.5cm	10	48
POREG85155	POREG85155	POREPUG85155	POREPUGFS85155	8.5cmx15.5cm	10	48
POREG85150	POREG85150	POREPUG85150	POREPUGFS85150	8.5cmx15cm	10	48
POREG8595	POREG8595	POREPUG8595	POREPUGFS8595	8.5cmx9.5cm	10	48
POREG90100	POREG90100	POREPUG90100	POREPUGFS90100	9cmx10cm	10	48
Protectfilm® Gentle	Protectfilm® GentleFrame Style		size	Pieces / box	Boxes / Case	
PROTG6070	PRGFS6070		6cmx7cm	100	16	
PROTG5075	PRGFS5075		5cmx7.5cm	100	16	
PROTG75100	PRGFS75100		7.5cmx10cm	100	16	
PROTG75120	PRGFS75120		7.5cmx12cm	100	16	
PROTG5082	PRGFS5082		5cmx8.2cm	100	16	
PROTG6585	PRGFS6585		6.5cmx8.5cm	100	16	
PROTG100120	PRGFS100120		10cmx12cm	50	16	
PROTG100140	PRGFS100140		10cmx14cm	50	16	
PROTG100160	PRGFS100160		10cmx16cm	50	16	
PROTG4545	PRGFS4545		4.5cmx4.5cm	100	16	
PROTG4867	PRGFS4867		4.8cmx6.7cm	100	16	
PROTG100100	PRGFS100100		10cmx10cm	100	16	
PROTG100150	PRGFS100150		10cmx15cm	50	16	
PROTG100200	PRGFS100200		10cmx20cm	40	16	
PROTG100250	PRGFS100250		10cmx25cm	40	16	
PROTG100300	PRGFS100300		10cmx30cm	40	16	
PROTG150200	PRGFS150200		15cmx20cm	50	8	
PROTG200300	PRGFS200300		20cmx30cm	25	8	
PROTG150150	PRGFS150150		15cmx15cm	25	8	
PROTG200200	PRGFS200200		20cmx20cm	25	8	
PROTG165180	PRGFS165180		16.5cmx18cm	25	8	
PROTG220220	PRGFS220220		22cmx22cm	25	8	
PROTG100180	PRGFS100180		10cmx18cm	50	16	
PROTG150250	PRGFS150250		15cmx25cm	25	8	
PROTG200250	PRGFS200250		20cmx25cm	25	8	
PROTG4870	PRGFS4870		4.8cmx7cm	100	16	
PROTG5070	PRGFS5070		5cmx7cm	100	16	
PROTG5072	PRGFS5072		5cmx7.2cm	100	16	
PROTG5080	PRGFS5080		5cmx8cm	100	16	
PROTG6080	PRGFS6080		6cmx8cm	100	16	
PROTG100350	PRGFS100350		10x35 cm			

Gentle Post-operative Dressing & OTC bandages

Sizes:

Cure-Aid® Gentle	size	Pieces / box	Boxes / Case
CUREG938	9 X 38 mm	100	60
CUREG1657	16 X 57 mm	100	60
CUREG1938	19 X 38 mm	100	60
CUREG1963	19 X 63 mm	100	60
CUREG1964	19 X 64 mm	100	60
CUREG1972	19 X 72 mm	100	60
CUREG1976	19 X 76 mm	100	60
CUREG2572	25 X 72 mm	100	60
CUREG2576	25 X 76 mm	100	60
CUREG3055	30 X 55 mm	100	60
CUREG3064	30 X 64 mm	50	60
CUREG3072	30 X 72 mm	50	60
CUREGANCHOR	30 X 72 mm Anchor	50	60
CUREG3838	38 X 38 mm	100	60
CUREG3872	38 X 72 mm	100	60
CUREG5072	50 X 72 mm	50	60
CUREG50100	50 X 100 mm	50	60
CUREG5060	50 X 60 mm Fingertip	50	60
CUREG4563	45 X 63 mm Butterfly	50	60
CUREG25	SPOT 25 mm	100	60
CUREG22	SPOT 22 mm	100	60
CUREG2525	25 X 25 mm	100	60
EYEPADG5062	50x62 mm Eyepad child	20	40
EYEPADG5882	58x82mm Eyepad Adult	20	40
EYEPADG6595	65x95 mm ocular	20	40
CUREGFEET	Feet	10	48
CUREG2235	22 x 35 mm	100	60
CUREG2665	26 x 65 mm	100	60
CUREG3070	30 x 70 mm	100	60
CUREG60100	60 x 100 mm	100	60
CUREG601000	60 x 1000 mm	100	60
CUREG64110	64 X 110 mm OVALADULT	100	60
CUREG5284	52 X 84 mm OVALCHILD	100	60
CUREG20120	20 x 120 mm	100	60
CUREG30120	30 x 120 mm	100	60
CUREG20180	20 x 180 mm	100	60
CUREG30180	30 x 180 mm	100	60
CUREG2582	25 x 82 mm	100	60
CUREG2956	29 x 56 mm	100	60
CUREGCHARACTERS	19x63 mm	100	60
CUREG2075	20x75 mm	100	60
CUREG1955	19 x 55 mm	100	60
CUREG4550	45 x 50 mm Butterfly	50	60
DERMIG60100	60 x 100 mm		
DERMIG601000	60 x 1000 mm		
DERMIG80100	80 X 100 mm		
DERMIG801000	80 x 1000 mm		
DERMIG40100	40 x 100 mm		
DERMIG60200	60 x 200 mm		
DERMIG2560	25 x 60 mm		
DERMIG401000	40 x 1000 mm		
DERMIG40500	40 x 500 mm		
DERMIG60500	60 x 500 mm		
DERMIG80500	80 x 500 mm		
DERMIG405000	40 x 5000 mm		
DERMIG605000	60 x 5000 mm		
DERMIG805000	80 x 5000 mm		
DERMIG20120	20 x 120 mm		
DERMIG30120	30 x 120 mm		
DERMIG20180	20 x 180 mm		

Gentle Post-operative Dressing & OTC bandages

Sizes:

Cure-Aid® Gentle	size	Pieces / box	Boxes / Case
DERMIG30180	30 x 180 mm		
DERMIG20120IW	20 x 120 mm individually wrapped		
DERMIG30120IW	30 x 120 mm individually wrapped		
DERMIG20180IW	20 x 180 mm individually wrapped		
DERMIG30180IW	30 x 180 mm individually wrapped		
DERMIG2040	20x40 mm		
DERMIG2060	20x60 mm		
DERMIG70100	70x100 mm		

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Gentle Post-Operative Dressing & OTC bandages						
Issue date : 22-02-2018		Prepared by : Marina Medhat						
Colors : CMYK		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	06-06-2022	Updating codes	File maintenance					