

Pharmacool Instant Ice Pack



Description

A portable, easy to use cold compress, that is highly effective in providing cold therapy for first aid and pain relief without the need of prior freezing or refrigerating.

Mode of Action

The instant ice pack is a fast and effective cooling method. It is easily activated by simply squeezing the ice pack, so it immediately becomes cold. This occurs due to bursting of the inner water bag which then becomes in contact with urea causing an endothermic reaction, so cooling occurs.

Features

- Instant activation with no freezing or refrigeration required.
- Instant cold therapy for first aid and pain relief.
- Convenient and portable size.
- Safe and cost effective.
- Flexible design.
- Easy to use.
- Disposable.
- Latex free.


Intended Purpose

It is used as first aid for sport injuries, to improve sprains, strains, bruises, muscle aches, decrease temperature by physiotherapy, tension headaches, minor burns. Also, it is used to reduce swelling and relieve discomfort from insect bites.

Instant Ice Pack	Size	Pieces/Box	Boxes/Case	Case Dimensions
ICE130150	13 x 15 cm	2	24	

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Instant Ice Pack						
Issue date : 21-02-2023		Prepared by : Nader Nabil Zaki						
Colors : 		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
1	N/A	New Issue	N/A			Artwork manager	Marketing department	RA Department
2	26-06-2023	Adding brand name	File maintenance					