

Medical Adhesive Tapes



Description:

Medical adhesive tapes are made of different backing materials, which have hypoallergenic adhesive to be soft on the skin and meet consumer requirements.

a. **Airopplast®**: is a water-resistant, transparent, perforated plastic medical tape that has a bi-directional perforation for easy tear by hand, coated on one side with hypoallergenic adhesive.

Also available non-perforated.

b. **Neoplast®**: is a water-repellent textile fabric (white/skin color), easy to tear off; coated on one side with pressure-sensitive adhesive permeable to air and water vapor, adheres reliably, and can be removed easily and leaves no residues on the skin.





Available in 100% cotton or 100% polyester.

c. **Pharmamed®**: is a microporous non-woven backing material coated on one side with a hypoallergenic adhesive.

d. **Pharmafix®**: is made of soft, breathable, spunlaced non-woven coated one side with a hypoallergenic adhesive and protected with grid printed silicon paper with a crack-back cut for easy application.

e. **Pharmafix® PU**: is made of thin, flexible, water-proof, bacterial & virus-proof polyurethane film coated on one side with a hypoallergenic adhesive and protected with a grid printed silicon paper with a crack-back cut for easy application.

f. **Silkplast®**: is a silk-like cloth tape made of either 100% cellulose acetate or 100% polyester coated on one side with pressure-sensitive adhesive.



Medical Adhesive Tapes

Intended Purpose:

Airoplast®:

- For dressing retention.
- Mainly used on patients who are sensitive to traditional adhesive tapes.

Pharmamed®:

- For skin sensitive to traditional adhesive tape.
- Babies with fragile skin.
- Elderly patients with dry skin.
- Patients undergoing radiation treatments.

Neoplast®:

- Provides preventive support for ligaments during sports.
- Water repellent.

Silkplast®:

- For dressings and medical device retention for short periods on skin.
- Fixation of needles, probes, catheters and medical equipment.
- Fixation of dressings and appliances where they are likely to be in contact with dirt or water.

Pharmafix®:

- Retention of post-operative wound dressings.
- Retaining gauze and absorbent pads.
- To secure catheters and drainage tubes.

Pharmafix® PU:

- Retention of primary dressing or low adherent dressings.
- To secure catheters and drainage tubes.
- Attachment of leg bags, stoma devices, etc.
- Breathable surgical tape.

Medical Adhesive Tapes

Sizes

Airoplast® Perforated	Airoplast® Non-Perforated	cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
AIRO124550	AIRO124550N	1.25cmx4.55m	½inx5y	24	12	
AIRO254550	AIRO254550N	2.5cmx4.55m	1inx5y	12	12	
AIRO504550	AIRO504550N	5cmx4.55m	2inx5y	6	12	
AIRO754550	AIRO754550N	7.5cmx4.55m	3inx5y	3	12	
AIRO1004550	AIRO1004550N	10cmx4.55m	4inx5y	2	12	
AIRO125000	AIRO125000N	1.25cmx5m	½inx5.47y	24	12	
AIRO255000	AIRO255000N	2.5cmx5m	1inx5.47y	12	12	
AIRO505000	AIRO505000N	5cmx5m	2inx5.47y	6	12	
AIRO755000	AIRO755000N	7.5cmx5m	3inx5.47y	3	12	
AIRO1005000	AIRO1005000N	10cmx5m	4inx5.47y	2	12	
AIRO129140	AIRO129140N	1.25cmx9.14m	½inx10y	24	12	
AIRO259140	AIRO259140N	2.5cmx9.14m	1inx10y	12	12	
AIRO509140	AIRO509140N	5cmx9.14m	2inx10y	6	12	
AIRO759140	AIRO759140N	7.5cmx9.14m	3inx10y	4	12	
AIRO1009140	AIRO1009140N	10cmx9.14m	4inx10y	3	12	
AIRO1210000	AIRO1210000N	1.25cmx10m	½inx10.94y	24	12	
AIRO2510000	AIRO2510000N	2.5cmx10m	1inx10.94y	12	12	
AIRO5010000	AIRO5010000N	5cmx10m	2inx10.94y	6	12	
AIRO7510000	AIRO7510000N	7.5cmx10m	3inx10.94y	4	12	
AIRO10010000	AIRO10010000N	10cmx10m	4inx10.94y	3	12	
AIROSET	AIROSETN	Assorted			12	
AIRO10015000		10cmX15m				
AIRO258000	AIRO258000N	2.5cmx8m				

Medical Adhesive Tapes

Silkplast® Latex	Silkplast® Latex Free	cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
SILK124550	SILKLF124550	1.25cmx4.55m	½inx5y	24	12	
SILK254550	SILKLF254550	2.5cmx4.55m	1inx5y	12	12	
SILK504550	SILKLF504550	5cmx4.55m	2inx5y	6	12	
SILK754550	SILKLF754550	7.5cmx4.55m	3inx5y	3	12	
SILK1004550	SILKLF1004550	10cmx4.55m	4inx5y	2	12	
SILK125000	SILKLF125000	1.25cmx5m	½inx5.47y	24	12	
SILK255000	SILKLF255000	2.5cmx5m	1inx5.47y	12	12	
SILK505000	SILKLF505000	5cmx5m	2inx5.47y	6	12	
SILK755000	SILKLF755000	7.5cmx5m	3inx5.47y	3	12	
SILK1005000	SILKLF1005000	10cmx5m	4inx5.47y	2	12	
SILK129140	SILKLF129140	1.25cmx9.14m	½inx10y	24	12	
SILK259140	SILKLF259140	2.5cmx9.14m	1inx10y	12	12	
SILK509140	SILKLF509140	5cmx9.14m	2inx10y	6	12	
SILK759140	SILKLF759140	7.5cmx9.14m	3inx10y	4	12	
SILK1009140	SILKLF1009140	10cmx9.14m	4inx10y	3	12	
SILK1210000	SILKLF1210000	1.25cmx10m	½inx10.94y	24	12	
SILK2510000	SILKLF2510000	2.5cmx10m	1inx10.94y	2	12	
SILK5010000	SILKLF5010000	5cmx10m	2inx10.94y	1	12	
SILK7510000	SILKLF7510000	7.5cmx10m	3inx10.94y	1	12	
SILK10010000	SILKLF10010000	10cmx10m	4inx10.94y	3	12	
SILKSET	SILKLFSET	Assorted			12	
SILK1805000	SILKLF1805000	18cmx5m		1	50	
SILK1005000P	SILKLF1005000P	10cmx5m Perforated	4inx5.47y	1	100	
SILK1805000P	SILKLF1805000P	18cmx5m Perforated	7inx4.37y	1	50	
SILK1004000P	SILKLF1004000P	10cmx4m Perforated	4inx5.47y	1	100	
SILK1804000P	SILKLF1804000P	18cmx4m Perforated	7inx4.37y	1	50	
Silkplast® PET Latex Free		cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
SILKLPET124550		1.25cmx4.55m	½inx5y	24	12	
SILKLPET254550		2.5cmx4.55m	1inx5y	12	12	
SILKLPET504550		5cmx4.55m	2inx5y	6	12	
SILKLPET754550		7.5cmx4.55m	3inx5y	3	12	
SILKLPET1004550		10cmx4.55m	4inx5y	2	12	
SILKLPET125000		1.25cmx5m	½inx5.47y	24	12	
SILKLPET255000		2.5cmx5m	1inx5.47y	12	12	
SILKLPET505000		5cmx5m	2inx5.47y	6	12	
SILKLPET755000		7.5cmx5m	3inx5.47y	3	12	
SILKLPET1005000		10cmx5m	4inx5.47y	2	12	
SILKLPET129140		1.25cmx9.14m	½inx10y	24	12	
SILKLPET259140		2.5cmx9.14m	1inx10y	12	12	
SILKLPET509140		5cmx9.14m	2inx10y	6	12	
SILKLPET759140		7.5cmx9.14m	3inx10y	4	12	

Medical Adhesive Tapes

Silkplast® PET Latex Free		cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
SILKLPET1009140		10cmx9.14m	4inx10y	3	12	
SILKLPET1210000		1.25cmx10m	½inx10.94y	24	12	
SILKLPET2510000		2.5cmx10m	1inx10.94y	2	12	
SILKLPET5010000		5cmx10m	2inx10.94y	1	12	
SILKLPET7510000		7.5cmx10m	3inx10.94y	1	12	
SILKLPET10010000		10cmx10m	4inx10.94y	3	12	
SILKLPETSET		Assorted			12	
SILKLPET1805000		18cmx5m		1	50	
SILKLPET1005000P		10cmx5m Perforated	4inx5.47y	1	100	
SILKLPET1805000P		18cmx5m Perforated	7inx4.37y	1	50	
SILKLPET1004000P		10cmx4m Perforated	4inx5.47y	1	100	
SILKLPET1804000P		18cmx4m Perforated	7inx4.37y	1	50	
Pharmafix®	Pharmafix® PU	cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
FIX254550	FIXPU254550	2.5cmx4.55m	1inx5y	2	72	
FIX504550	FIXPU504550	5cmx4.55m	2inx5y	1	72	
FIX754550	FIXPU754550	7.5cmx4.55m	3inx5y	1	36	
FIX1004550	FIXPU1004550	10cmx4.55m	4inx5y	1	36	
FIX1504550	FIXPU1504550	15cmx4.55m	6inx5y	1	36	
FIX2004550	FIXPU2004550	20cmx4.55m	8inx5y	1	18	
FIX40150	FIXPU40150	4 cm x 15 cm	1.6inx6in	100	36	

Medical Adhesive Tapes

Pharmafix®	Pharmafix® PU	cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
FIX3004550	FIXPU3004550	30cmx4.55m	12inx5y	1	18	
FIX255000	FIXPU255000	2.5cmx5m	1inx5.47y	2	72	
FIX505000	FIXPU505000	5cmx5m	2inx5.47y	1	72	
FIX755000	FIXPU755000	7.5cmx5m	3inx5.47y	1	36	
FIX1005000	FIXPU1005000	10cmx5m	4inx5.47y	1	36	
FIX1505000	FIXPU1505000	15cmx5m	6inx5.47y	1	36	
FIX2005000	FIXPU2005000	20cmx5m	8inx5.47y	1	18	
FIX3005000	FIXPU3005000	30cmx5m	12inx5.47y	1	18	
FIX259140	FIXPU259140	2.5cmx9.14m	1inx10y	2	72	
FIX509140	FIXPU509140	5cmx9.14m	2inx10y	1	72	
FIX759140	FIXPU759140	7.5cmx9.14m	3inx10y	1	36	
FIX1009140	FIXPU1009140	10cmx9.14m	4inx10y	1	36	
FIX1509140	FIXPU1509140	15cmx9.14m	6inx10y	1	36	
FIX2009140	FIXPU2009140	20cmx9.14m	8inx10y	1	18	
FIX3009140	FIXPU3009140	30cmx9.14m	12inx10y	1	18	
FIX2510000	FIXPU2510000	2.5cmx10m	1inx10.94y	2	72	
FIX5010000	FIXPU5010000	5cmx10m	2inx10.94y	1	72	
FIX7510000	FIXPU7510000	7.5cmx10m	3inx10.94y	1	36	
FIX10010000	FIXPU10010000	10cmx10m	4inx10.94y	1	36	
FIX15010000	FIXPU15010000	15cmx10m	6inx10.94y	1	36	
FIX20010000	FIXPU20010000	20cmx10m	8inx10.94y	1	18	
FIX30010000	FIXPU30010000	30cmx10m	12inx10.94y	1	18	
FIX1805000		18cmx5m	7inx5.47y	1	18	
FIXSET	FIXPUSET	Assorted			36	
FIX20150		2 cmx 15 cm	0.8in x 6in	100	36	
FIX1801000		18cmx1m	7in x 1y	1	36	
Pharmamed®		cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
MED124550		1.25cmx4.55m	½inx5y	24	12	
MED254550		2.5cmx4.55m	1inx5y	12	12	
MED504550		5cmx4.55m	2inx5y	6	12	
MED754550		7.5cmx4.55m	3inx5y	3	12	
MED1004550		10cmx4.55m	4inx5y	2	12	
MED125000		1.25cmx5m	½inx5.47y	24	12	
MED255000		2.5cmx5m	1inx5.47y	12	12	
MED505000		5cmx5m	2inx5.47y	6	12	
MED755000		7.5cmx5m	3inx5.47y	3	12	
MED1005000		10cmx5m	4inx5.47y	2	12	
MED129140		1.25cmx9.14m	½inx10y	24	12	
MED259140		2.5cmx9.14m	1inx10y	12	12	
MED509140		5cmx9.14m	2inx10y	6	12	
MED759140		7.5cmx9.14m	3inx10y	4	12	

Medical Adhesive Tapes

Neoplast® PET Latex free	cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
NEOLFPET124550	1.25cmx4.55m	½inx5y	24	12	
NEOLFPET254550	2.5cmx4.55m	1inx5y	12	12	
NEOLFPET504550	5cmx4.55m	2inx5y	6	12	
NEOLFPET754550	7.5cmx4.55m	3inx5y	3	12	
NEOLFPET1004550	10cmx4.55m	4inx5y	2	12	
NEOLFPET125000	1.25cmx5m	½inx5.47y	24	12	
NEOLFPET255000	2.5cmx5m	1inx5.47y	12	12	
NEOLFPET505000	5cmx5m	2inx5.47y	6	12	
NEOLFPET755000	7.5cmx5m	3inx5.47y	3	12	
NEOLFPET1005000	10cmx5m	4inx5.47y	2	12	
NEOLFPET129140	1.25cmx9.14m	½inx10y	24	12	
NEOLFPET259140	2.5cmx9.14m	1inx10y	12	12	
NEOLFPET509140	5cmx9.14m	2inx10y	6	12	
NEOLFPET759140	7.5cmx9.14m	3inx10y	4	12	
NEOLFPET1009140	10cmx9.14m	4inx10y	3	12	
NEOLFPET1210000	1.25cmx10m	½inx10.94y	24	12	
NEOLFPET2510000	2.5cmx10m	1inx10.94y	12	12	
NEOLFPET5010000	5cmx10m	2inx10.94y	6	12	
NEOLFPET7510000	7.5cmx10m	3inx10.94y	4	12	
NEOLFPET10010000	10cmx10m	4inx10.94y	3	12	
NEOLFPETSET	Assorted			12	
NEOLFPET1805000	18cmx5m	7inx4.37y	1	50	
NEOLFPET1005000P	10cmx5m Perforated	4inx5.47y	1	100	
NEOLFPET1805000P	18cmx5m Perforated	7inx4.37y	1	50	
NEOLFPET1004000P	10cmx4m Perforated	4inx5.47y	1	100	
NEOLFPET1804000P	18cmx4m Perforated	7inx4.37y	1	50	
NEOLFPET3710000	3.75cm x 10m	1.4in x 10.94 Y	6	12	

Medical Adhesive Tapes

Pharmamed®		cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
MED1009140		10cmx9.14m	4inx10y	3	12	
MED1210000		1.25cmx10m	½inx10.94y	4	72	
MED2510000		2.5cmx10m		2	72	
MED5010000		5cmx10m	1inx10.94y	1	72	
MED7510000		7.5cmx10m	2inx10.94y	1	36	
MED10010000		10cmx10m	3inx10.94y	1	36	
MEDSET		Assorted	4inx10.94y			
Neoplast® Latex	Neoplast® Latex free	cm	Inch	Rolls / box	Boxes / Case	Case Dimensions
NEO124550	NEOLF124550	1.25cmx4.55m	½inx5y	24	12	
NEO254550	NEOLF254550	2.5cmx4.55m	1inx5y	12	12	
NEO504550	NEOLF504550	5cmx4.55m	2inx5y	6	12	
NEO754550	NEOLF754550	7.5cmx4.55m	3inx5y	3	12	
NEO1004550	NEOLF1004550	10cmx4.55m	4inx5y	2	12	
NEO125000	NEOLF125000	1.25cmx5m	½inx5.47y	24	12	
NEO255000	NEOLF255000	2.5cmx5m	1inx5.47y	12	12	
NEO505000	NEOLF505000	5cmx5m	2inx5.47y	6	12	
NEO705000	NEOLF705000	7cmx5m	2.8inx5.47y	3	12	
NEO755000	NEOLF755000	7.5cmx5m	3inx5.47y	3	12	
NEO1005000	NEOLF1005000	10cmx5m	4inx5.47y	2	12	
NEO129140	NEOLF129140	1.25cmx9.14m	½inx10y	24	12	
NEO259140	NEOLF259140	2.5cmx9.14m	1inx10y	12	12	
NEO509140	NEOLF509140	5cmx9.14m	2inx10y	6	12	
NEO759140	NEOLF759140	7.5cmx9.14m	3inx10y	4	12	
NEO1009140	NEOLF1009140	10cmx9.14m	4inx10y	3	12	
NEO1210000	NEOLF1210000	1.25cmx10m	½inx10.94y	24	12	
NEO2510000	NEOLF2510000	2.5cmx10m	1inx10.94y	12	12	
NEO5010000	NEOLF5010000	5cmx10m	2inx10.94y	6	12	
NEO7510000	NEOLF7510000	7.5cmx10m	3inx10.94y	4	12	
NEO10010000	NEOLF10010000	10cmx10m	4inx10.94y	3	12	
NEOSET	NEOLFSET	Assorted			12	
NEO1805000	NEOLF1805000	18cmx5m	7inx4.37y	1	50	
NEO1005000P	NEOLF1005000P	10cmx5m Perforated	4inx5.47y	1	100	
NEO1805000P	NEOLF1805000P	18cmx5m Perforated	7inx4.37y	1	50	
NEO1004000P	NEOLF1004000P	10cmx4m Perforated	4inx5.47y	1	100	
NEO1804000P	NEOLF1804000P	18cmx4m Perforated	7inx4.37y	1	50	
NEO3710000	NEOLF3710000	3.75cm x 10m	1.4in x 10.94 Y	6	12	
NEO3712000	NEOLF3712000	3.75cm x 12m	1.4in x 13.12y	6	12	

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Medical Adhesive Tapes						
Issue date : 02-04-2020		Prepared by : Diana Amr						
Colors : <input type="checkbox"/>		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	31-07-2019	- changing Pharmafix PU pic via request from QA - Waterproof Pharmamed removed	CRS Request					
3	16-09-2019	Adding extra sizes for table of sizes	Marketing point of view					
4	02-04-2020	- adding case dimensions - modifying some numbers of units - adding new sizes	- file maintenance - file update					
5	23-6-2021	adding new codes						
6	27-09-2021	aminding codes	- file maintenance - file update					
7	06-02-2022	aminding codes	- file maintenance					
8	02-02-2023	adding new codes adding silkplast PET latex variant	- file maintenance					
9	21-05-2023	Text amendments	- file maintenance					
10	28-05-2023	Text update Adding new size	- file maintenance					
11	18-03-2024	Adding new size	- file maintenance					
12	01-04-2024	Adding new size	- file maintenance					
13	16-04-2024	all waterproof are removed	- file maintenance					
14	19-05-2024	Adding new size	- file maintenance					
15	09-10-2024	adding neoplast PET Latex Free Silkplast Pet latex removed	- file maintenance					