

# Pharma – Undercast Pad <sup>®</sup>



## Description

- Synthetic Orthopaedic Padding

## Features

- Excellent conformability
- Easy to tear
- Autoclavable
- Good cushioning
- Soft
- Polyester and acrylic fibre mix

## Benefits

- Due to excellent conformability allows even coverage and a smooth finish

## Indications

- Suitable under traditional plaster of Paris casts and all types of synthetic casts in the treatment of fractures or orthopaedic conditions
- Also suitable for use with splinting material or other situations where padding is required
- Provides cushioned protection of the injured limb

Product code	Size	Rolls / case	Case dimensions
CAST502700	5x270cm	200	57 x 32 x 28 cm
CAST752700	7.5x270cm	150	57 x 32 x 28 cm
CAST1002700	10x270cm	100	57 x 32 x 28 cm
CAST1502700	15x270cm	75	57 x 32 x 28 cm
CAST2002700	20x270cm	50	57 x 32 x 28 cm

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263

- Kafr-El-Zayat P.O.Box 38 Egypt

Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail :sales@pharmaplast-online.com



www.pharmaplast-online.com

Pharmaplast		Amria 23612, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3 )				
Customer : Pharmaplast				Product name : Pharma – Undercast Pad <sup>®</sup>				
Issue date : 19-07-2015				Prepared by : Nader Nabil Zaki				
Colors :				Dimension : A4				
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	16-05-2022	Code updaed	File maintenance					