

# Pharmacoll® Family



## Description:

The **Pharmacoll®** product line represents a gelatine free hydrocolloid wound dressing containing Carboxymethyl Cellulose powder that converts into a gel when it comes in contact with wound fluid. The moist wound environment created under the dressing speeds the natural healing process, reduces pain, and protects against bacterial and other external contaminants.

- **Pharmacoll®, Pharmacoll Thin and Pharmacoll® Extra Thin** are gelatin free hydrocolloid adhesive wound dressings, laminated to a semi-permeable, bacteria-proof, virus-proof and waterproof polyurethane film backing. The difference between the 3 variants is the thickness of the hydrocolloid layer, where **Pharmacoll®** is 1 mm thickness, **Pharmacoll® Thin** 0.6 mm and **Pharmacoll® Extra Thin** is 0.3 mm.
- **Pharmacoll® Cushion** is a gelatin free hydrocolloid adhesive wound dressing laminated to a waterproof thick cross-linked polyethylene foam film backing.
- **Pharmacoll® Comfort** is a gelatin free hydrocolloid adhesive wound dressing laminated to a semi-permeable, bacteria-proof, virus-proof and waterproof polyurethane film backing with adhesive polyurethane film border.
- **Pharmacoll® Comfort Plus** is thick hydrocolloid island pad placed in the middle of a thin hydrocolloid layer. The thin hydrocolloid layer is laminated with a semi permeable, bacteria proof, viral proof and waterproof polyurethane film.

## Mode of Action:

When **Pharmacoll®** dressing comes into contact with wound exudates, the hydrocolloid swells and forms a gel that covers the wound and maintains moisture. This gel continues absorbing the wound fluid and minimizes adherence to the wound. When the dressing is changed, the wound stays protected by a film of gel.

## Features:

- Adhere to both dry and moist skin surfaces.
- Deliver gentle and secure adhesion that does not irritate the skin.
- Have fluid absorption capacity that is adjustable based on adhesive formula.
- Provide barrier properties to external contaminants with secure adhesion to skin.
- Are flexible and conformable and mold to difficult anatomical contours.
- Cushion effect.
- Latex free.
- Tapered edge minimizes bunching and rolling up extending wear time and lowering treatment cost.
- A tapered hydrocolloid dressing is thin at the edges and thicker in the center. Because of its thin edges, tape is not required and the dressing is more flexible and conformable.

## Intended Purpose:

**Pharmacoll®** dressings are indicated to be used in:

- Partial- and full-thickness dermal ulcers.
- Pressure, Leg and Diabetic ulcers.
- Burns (First and Second Degrees).
- Skin Donor sites and Skin grafts.
- Superficial wounds/minor abrasions.
- Protection of fragile skin or areas exposed to friction such as elbows, heels or the sacral area.
- Ostomy wounds to protect the skin from stomal effluent.
- Post-op dressings for suture line protection.

# Pharmacoll® Family

Pharmacoll®	Size	Pieces/ box	Boxes / Case	Case Dimensions
COLL5050	50x50mm	10	48	37x21x16
COLL5075	50x75mm	10	48	37x21x16
COLL50100	50x100mm	10	48	37x21x16
COLL7575	75x75mm	10	48	37x21x16
COLL100100	100x100mm	10	48	37x21x16
COLL125125	125x125mm	5	24	38x22x25
COLL100200	100x200mm	5	24	52x30x26
COLL150150	150x150mm	5	24	38x22x25
COLL150200	150x200mm	5	24	52x30x26
COLL200200	200x200mm	5	24	52x30x26
COLL200300	200x300mm	5	24	52x30x36
COLL75100	75x100mm	10	48	37x21x16
COLL67110	67x110mm	10	48	37x21x16
COLL4368	43x68mm	10	48	37x21x16
COLL2868	28x68mm	10	48	37x21x16
COLL4565	45x65mm	10	48	37x21x16
COLL165180	165x180mm concave	5	24	44x25x29
COLL220220	220x220mm sacral	5	24	64x36x32
COLL170170	170x170mm sacral	5	24	44x25x29
COLL180180	180x180mm	5	24	52x30x26
COLL50150	50x150mm	5	24	38x22x25
COLL50200	50x200mm	5	24	52x30x26
COLL50250	50x250mm	5	24	52x30x36
COLL2045	20x45mm	10	48	37x21x16
COLL5072	50x72mm	10	48	37x21x16
COLL100150	100x150mm	5	24	52x30x26
COLL3844	38X44 mm	10	48	37x21x16
COLL100120	100x120 mm	5	24	38x22x25
COLL130150	130x150 mm	5	24	38x22x25
COLL6080	60x80 mm	10	48	37x21x16
COLL1501000	150x1000 mm	5	24	
COLL100350	100x350 mm	5	24	
COLL100250	100x250 mm	5	24	

# Pharmacoll® Family

Pharmacoll® Thin	Size	Pieces/ box	Boxes / Case	Case Dimensions
COLLT5050	50x50mm	10	48	37x21x16
COLLT1501000	150 mm x 1 mt	1	48	
COLLT5075	50x75mm	10	48	37x21x16
COLLT7575	75x75mm	10	48	37x21x16
COLLT50100	50x100mm	10	48	37x21x16
COLLT100100	100x100mm	10	48	37x21x16
COLLT125125	125x125mm	5	24	38x22x25
COLLT100200	100x200mm	5	24	52x30x26
COLLT150150	150x150mm	5	24	38x22x25
COLLT150200	150x200mm	5	24	52x30x26
COLLT200200	200x200mm	5	24	52x30x26
COLLT200300	200x300mm	5	24	52x30x36
COLLT75100	75x100mm	10	48	37x21x16
COLLT67110	67x110mm	10	48	37x21x16
COLLT4368	43x68mm	10	48	37x21x16
COLLT2868	28x68mm	10	48	37x21x16
COLLT4565	45x65mm	10	48	37x21x16
COLLT165180	165x180mm concave	5	24	44x25x29
COLLT220220	220x220mm sacral	5	24	64x36x32
COLLT170170	170x170mm sacral	5	24	44x25x29
COLLT180180	180x180mm	5	24	52x30x26
COLLT50150	50x150mm	5	24	38x22x25
COLLT50200	50x200mm	5	24	52x30x26
COLLT50250	50x250mm	5	24	52x30x36
COLLT100150	100x150mm	5	24	38x22x25
COLLT2045	20x45mm	10	48	37x21x16
COLLT5072	50x72mm	10	48	37x21x16
COLLT3844	38X44 mm	10	48	37x21x16
COLLT100120	100x120 mm	5	24	38x22x25
COLLT130150	130x150 mm	5	24	38x22x25
COLLT100350	100x350 mm	5	24	
COLLT100250	100x250 mm	5	24	
COLLT6080	60x80 mm	10	48	

# Pharmacoll® Family

Pharmacoll® ExtraThin	Size	Pieces/ box	Boxes / Case	Case Dimensions
COLLEXTRAT5050	50x50mm	10	48	37x21x16
COLLEXTRAT5075	50x75mm	10	48	37x21x16
COLLEXTRAT7575	75x75mm	10	48	37x21x16
COLLEXTRAT50100	50x100mm	10	48	37x21x16
COLLEXTRAT100100	100x100mm	10	48	37x21x16
COLLEXTRAT125125	125x125mm	5	24	38x22x25
COLLEXTRAT100200	100x200mm	5	24	52x30x26
COLLEXTRAT150150	150x150mm	5	24	38x22x25
COLLEXTRAT150200	150x200mm	5	24	52x30x26
COLLEXTRAT200200	200x200mm	5	24	52x30x26
COLLEXTRAT200300	200x300mm	5	24	52x30x36
COLLEXTRAT75100	75x100mm	10	48	37x21x16
COLLEXTRAT67110	67x110mm	10	48	37x21x16
COLLEXTRAT4368	43x68mm	10	48	37x21x16
COLLEXTRAT2868	28x68mm	10	48	37x21x16
COLLEXTRAT4565	45x65mm	10	48	37x21x16
COLLEXTRAT165180	165x180mm concave	5	24	44x25x29
COLLEXTRAT220220	220x220mm sacral	5	24	64x36x32
COLLEXTRAT170170	170x170mm sacral	5	24	44x25x29
COLLEXTRAT180180	180x180mm	5	24	52x30x26
COLLEXTRAT50150	50x150mm	5	24	38x22x25
COLLEXTRAT50200	50x200mm	5	24	52x30x26
COLLEXTRAT50250	50x250mm	5	24	52x30x36
COLLEXTRAT100150	100x150mm	5	24	38x22x25
COLLEXTRAT2045	20x45mm	10	48	37x21x16
COLLEXTRAT5072	50x72mm	10	48	37x21x16
COLLEXTRAT3844	38X44 mm	10	48	37x21x16
COLLEXTRAT100120	100x120 mm	5	24	38x22x25
COLLEXTRAT130150	130x150 mm	5	24	38x22x25
COLLEXTRAT6080	60x80 mm	10	48	
COLLEXTRAT1501000	150x1000 mm	5	24	
COLLEXTRAT100350	100x350 mm	5	24	
COLLEXTRAT100250	100x250 mm	5	24	

# Pharmacoll® Family

Pharmacoll® Cushion	Size	Pieces/ box	Boxes / Case	Case Dimensions
COLLPEF5050	50x50mm	10	48	37x21x16
COLLPEF5075	50x75mm	10	48	37x21x16
COLLPEF7575	75x75mm	10	48	37x21x16
COLLPEF50100	50x100mm	10	48	37x21x16
COLLPEF100100	100x100mm	10	48	37x21x16
COLLPEF125125	125x125mm	5	24	38x22x25
COLLPEF100200	100x200mm	5	24	52x30x26
COLLPEF150150	150x150mm	5	24	38x22x25
COLLPEF150200	150x200mm	5	24	52x30x26
COLLPEF200200	200x200mm	5	24	52x30x26
COLLPEF200300	200x300mm	5	24	52x30x36
COLLPEF75100	75x100mm	10	48	37x21x16
COLLPEF67110	67x110mm	10	48	37x21x16
COLLPEF4368	43x68mm	10	48	37x21x16
COLLPEF2868	28x68mm	10	48	37x21x16
COLLPEF4565	45x65mm	10	48	37x21x16
COLLPEF165180	165x180mm concave	5	24	44x25x29
COLLPEF220220	220x220mm sacral	5	24	64x36x32
COLLPEF170170	170x170mm sacral	5	24	44x25x29
COLLPEF180180	180x180mm	5	24	52x30x26
COLLPEF50150	50x150mm	5	24	38x22x25
COLLPEF50200	50x200mm	5	24	52x30x26
COLLPEF50250	50x250mm	5	24	52x30x36
COLLPEF100150	100x150mm	5	24	38x22x25
COLLPEF2045	20x45mm	10	48	37x21x16
COLLPEF5072	50x72mm	10	48	37x21x16
COLLPEF3844	38X44 mm	10	48	37x21x16
COLLPEF100120	100x120 mm	5	24	38x22x25
COLLPEF130150	130x150 mm	5	24	38x22x25
COLLPEF6080	60x80 mm	10	48	
COLLPEF1501000	150x1000 mm	5	24	
COLLPEF100350	100x350 mm	5	24	
COLLPEF100250	100x250 mm	5	24	

# Pharmacoll® Family

Pharmacoll® comfort	Size	Pieces/ box	Boxes / Case	Case Dimensions
COLLC7575	75x75mm	10	48	37x21x16
COLLC75100	75x100mm	10	48	37x21x16
COLLC100100	100x100mm	10	48	37x21x16
COLLC200200	200x200mm	5	24	52x30x26
COLLC125125	125x125mm	5	24	38x22x25
COLLC150150	150x150mm	5	24	38x22x25
COLLC180180	180x180mm	5	24	52x30x26
COLLC4565	45x65mm	10	48	37x21x16
COLLC/PET4368	43x68mm	6	48	37x21x16
COLLC4368	43x68mm	6		37x21x16
		10		37x21x16
COLLC2045	20x45mm	10	48	37x21x16
COLLC/PET2045	20x45mm	6		37x21x16
COLLC165180	165x180mm concave	5	24	44x25x29
COLLC220220	220x220mm sacral	5	24	64x36x32
COLLC170170	170x170mm sacral	5	24	44x25x29
COLLC2235	22x35mm	10	48	37x21x16
COLLC2665	26x65mm	10	48	37x21x16
COLLC3070	30x70mm	10	48	37x21x16
COLLC60100	60x100mm	10	48	37x21x16
COLLC100120	100x120 mm	5	24	38x22x25
COLLC100130	100x130 mm	5	24	38x22x25
COLLC130150	130x150 mm	5	24	38x22x25
COLLC160170	160x170 mm	5	24	
COLLC67110	67x110 mm	10	48	
COLLC5050	50x50 mm	10	48	

# Pharmacoll® Family

Pharmacoll® Comfort Plus	size	Pieces/box	Boxes / case	Case Dimensions
COLLCPLUS7575	75x75mm	10	48	37x21x16
COLLCPLUS100100	100x100mm	10	48	37x21x16
COLLCPLUS125125	125x125mm	5	24	38x22x25
COLLCPLUS150150	150x150mm	5	24	38x22x25
COLLCPLUS180180	180x180mm	5	24	52x30x26
COLLCPLUS4565	45x65mm	10	48	37x21x16
COLLCPLUS4368	43x68mm	6	48	37x21x16
		10		37x21x16
COLLCPLUS2045	20x45mm	10	48	37x21x16
COLLCPLUS165180	165x180mm concave	5	24	44x25x29
COLLCPLUS220220	220x220mm sacral	5	24	64x36x32
COLLCPLUS170170	170x170mm sacral	5	24	44x25x29
COLLCPLUS2235	22x35mm	10	48	37x21x16
COLLCPLUS2665	26x65mm	10	48	37x21x16
COLLCPLUS3070	30x70mm	10	48	37x21x16
COLLCPLUS60100	60x100mm	10	48	37x21x16
COLLCPLUS67110	67x110mm	10	48	37x21x16
COLLCPLUS75100	75x100mm	10	48	37x21x16
COLLCPLUS200200	200x200mm	5	24	52x30x26
COLLCPLUS5050	50x50 mm	10	48	37x21x16
COLLCPLUS100120	100x120 mm	5	24	38x22x25
COLLCPLUS100130	100x130 mm	5	24	38x22x25
COLLCPLUS130150	130x150 mm	5	24	38x22x25
COLLCPLUS50100	50x100 mm	10	48	
COLLCPLUS160170	160x170 mm	5	24	

**- Amria - Alexandria - Egypt**

Tel. : 002 03 4500264 Fax : 002 03 4500263

**- Kafr-El-Zayat P.O.Box 38 Egypt**

Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com



www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Pharmacoll Family						
Issue date : 07/04/2020		Prepared by : Marina Medhat						
Colors : <input type="checkbox"/>		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	20-12-2021	Adding New Codes Adding case Dimensions						
3	11-03-2023	Updating t Adding New sizes						
4	09-04-2024	Text amendments	file maintenance					