


# Silicone adhesive tapes

Atraumatic fixation tapes with silicone adhesive.



## Description:

**a-Silotac tapes:** Waterproof rigid fabric coated with soft silicone adhesive and covered with release liner.

**b-Silofix tapes:** paper or non-woven backing material coated with soft silicone adhesive, for patients with sensitive skin.

## Mode of Action:


It has been shown that soft silicone adhesive does not cause trauma when removed from the skin. It has been described as 'atraumatic' for this reason.

Silicone adhesive removes less stratum corneum when compared to a composite hydrocolloid and polyurethane foam using an acrylic adhesive.


Soft Silicone adhesive is a material that adheres readily to intact dry skin but does not stick to the surface of a moist wound and does not cause damage upon removal.

## Features:

- Gentle Silicone adhesive.
- Suitable for fragile, sensitive skin or people allergic to traditional tapes or adhesives.
- Allows atraumatic dressing removal without leaving adhesive residue.
- May be lifted and reapplied if needed.
- Breathable and moisture-proof outermost film.
- Conformable.
- Unlike Silofix, Silotac is strong due to its rigid backing material.



No Skin stripping with Silotac



Painful skin stripping with traditional adhesive

## Intended Purpose:

It is used to secure wound dressings, catheters and other medical devices, gastrostomy tubes and other feeding tubes.

Silotac	Size	Units/ Box	Boxes/ Case
SILOT203000	2 x 300 cm	1	48
SILOT401500	4 x 150 cm	1	48
SILOT201500	2 x 150 cm	1	48
SILOT403000	4 x 300 cm	1	48
Silofix	Size	Units/ Box	Boxes/ Case
SILOF201500	2 x 150 cm	1	48
SILOF203000	2 x 300 cm	1	48
SILOF401500	4 x 150 cm	1	48
SILOF403000	4 x 300 cm	1	48
SILOF1254500	1.25 x 450 cm	1	48
SILOF254500	2.5 x 450 cm	1	48
SILOF504500	5 x 450 cm	1	48
SILOF754500	7.5 x 450 cm	1	48
SILOF1004500	10 x 450 cm	1	48
SILOF1255000	1.25 x 500 cm	1	48
SILOF255000	2.5 x 500 cm	1	48
SILOF505000	5 x 500 cm	1	48
SILOF755000	7.5 x 500 cm	1	48
SILOF1005000	10 x 500 cm	1	48

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263


- Kafr-El-Zayat P.O.Box 38 Egypt

Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com

 **Pharmaplast**

www.pharmaplast-online.com

Pharmaplast		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Silicone Tape						
Issue date : 31-01-2021		Prepared by : Nader Nabil Zaki						
Colors : 		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	28-07-2022	new codes	File maintenance					