

# PRODUCT BENEFITS

- Tack free: enhances unrolling tension and facilitates application and moulding with standard examination gloves.
- Porous: allows the skin to breathe.
- Radiolucent.
- Comfortable: Lightweight for easy mobility with a very smooth surface to enhance patient comfort.
- Conformable: stretches to conform to all body contours for optimal fit.
- Durable rigid support: long-lasting, high strength, yet lightweight cast, with high lamination between cast layers to improve cast strength.
- Offered in different colors to suit different patients' tastes.

## AVAILABLE SIZES

Sizes in cm	Sizes in inches
2.5 cm * 1.8 m	1 inch * 2 yards
5 cm * 3.6 m	2 inches * 4 yards
7.5 cm * 3.6 m	3 inches * 4 yards
10 cm * 3.6 m	4 inches * 4yards
12.5 cm * 3.6 m	5 inches * 4 yards

## INTENDED USE

Orthopedic casts for external rigid immobilization.



# SYNOCAST

Synocast is a lightweight durable fiberglass casting tape coated with a polyurethane resin.

When exposed to moisture or water, it hardens to provide excellent rigid support.



# SYNOCAST

## Fiberglass Cast