



VENDAPRESS FOAM

Vendapress foam is an elastic self-adherent foam bandage, that doesn't stick to skin nor hair, easy to apply without safety pins or fasteners.

Vendapress foam is made of soft conformable foam, that protects your wound from dirt and contaminants.

Available Colors



Blue



Beige

Contact




Sales@pharmaplast-online.net



www.pharmaplast-online.net



**Amria free zone 23512.
Alexandria Egypt.**



VENDAPRESS FOAM

Elastic Self-adherent Foam Bandage



Product Benefits

- Self-adherent: Sticks only to itself, does not stick to skin, wounds or hair, making bandage removal easy and painless.
- Easy to tear or cut.
- Soft, comfortable to wear and easy to use.
- Porous, allowing wound and skin to breathe.
- Water resistant: stays in place even when it gets wet.
- Provides moderate to high pressure to stop bleedings.
- Highly conformable, flexible and elastic: ideal for bandaging challenging body contours, e.g. finger, hand or toe.
- Available in roll form, it can be cut to fit different types of wounds, regardless of size.
- Available in 2 colors: Beige and blue:
- Beige: recommended for clean environment, easy to detect if bandage gets dirty or blood stained.
- Blue: recommended for food industry, highly visible and easy to detect.
- Available either with latex or latex free.
- Non-sterile.



AVAILABLE SIZES

3 * 450 cm
6 * 40 cm
6 * 100 cm
6 * 200 cm
6 * 450 cm
6 * 500 cm

Indented Use

- To protect minor wounds.
- Extra support or securement of padded plaster and dressings.

