


# eXsorb® and eXsorb® Carbon Family



## Description:

**eXsorb®** is a superabsorbent wound dressing enveloped in a non-woven pouch. Its white side is a low-adherent wound contact layer. Upon absorption of exudates, the dressing forms a soft hydrophilic gel that locks away bacteria and harmful proteases, reducing risk of infection and maceration.

**eXsorb® Algi** used with wounds that are prone to minor bleeding due to the haemostatic properties of the added calcium alginate.

**eXsorb® Carbon** same as **eXsorb®** but with an extra carbonized non-woven felt that is believed to adsorb offensive wound odour.

**eXsorb® island dressings** consist of superabsorbent pad (with or without carbon felt) laminated to a low-adherent wound contact layer and having an adhesive border. Upon absorption of exudates, the dressing forms a soft hydrophilic gel that locks away bacteria and harmful proteases, reducing risk of infection and maceration.

**eXsorb® OTC bandages** consist of superabsorbent pad with different backing materials like polyethylene film, polyurethane film, spunlaced nonwoven...etc.

## Intended Use:

Primary dressing indicated for partial and full thickness wounds, with moderate to heavy exudates that have breached skin dermis.

- **eXsorb® and eXsorb® island dressings:** Leg ulcers, pressure ulcers, diabetic ulcers, second degree surface burns, post-operative wounds.
- **eXsorb® Algi:** wounds that are prone to minor bleeding.
- **eXsorb® Carbon:** management of malodorous acute or chronic wounds.
- **eXsorb® OTC bandages** are used for minor wounds, cuts, grazes, minor abrasion.

## Features:

1. **High Absorption capacity:** Where it absorbs wound exudates and locks it inside the formed gel so prevents skin maceration.
2. **Conformable and Soft:** It conforms to wound contours due to Gel formation which prevents adherence of the dressing to the newly formed tissues
3. **Fewer dressing changes:** Minimizing treatment costs and disturbance to the wound.
4. **Reduction of Bioburden:** as it locks away bacteria and harmful proteases which promotes healing.
5. Suitable for use under compression.

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

eXsorb <sup>®</sup> Pad	eXsorb <sup>®</sup> Pad Carbon	eXsorb Algi	Size	Pieces/ box	Boxes/ Case	Case Dimensions
PADEX5050	PADEXC5050	PADEXALGI5050	5X5 cm	10	48	
PADEX5075	PADEXC5075	PADEXALGI5075	5X7.5 cm	10	24	
PADEX75100	PADEXC75100	PADEXALGI75100	7.5X10 cm	10	24	
PADEX75150	PADEXC75150	PADEXALGI75150	7.5X15 cm	5	24	
PADEX75200	PADEXC75200	PADEXALGI75200	7.5X20 cm	5	24	
PADEX50100	PADEXC50100	PADEXALGI50100	5X10 cm	10	24	
PADEX100100	PADEXC100100	PADEXALGI100100	10X10 cm	10	24	
PADEX100150	PADEXC100150	PADEXALGI100150	10X15 cm	5	24	
PADEX100200	PADEXC100200	PADEXALGI100200	10X20 cm	5	24	
PADEX100250	PADEXC100250	PADEXALGI100250	10X25 cm	5	24	
PADEX100300	PADEXC100300	PADEXALGI100300	10X30 cm	5	24	
PADEX1001000	PADEXC1001000	PADEXALGI1001000	10X100 cm	1	48	
PADEX1002000	PADEXC1002000	PADEXALGI1002000	10X200 cm	1	48	
PADEX150150	PADEXC150150	PADEXALGI150150	15X15 cm	5	24	
PADEX150200	PADEXC150200	PADEXALGI150200	15X20 cm	5	24	
PADEX150300	PADEXC150300	PADEXALGI150300	15X30 cm	5	24	
PADEX200200	PADEXC200200	PADEXALGI200200	20X20 cm	5	24	
PADEX200250	PADEXC200250	PADEXALGI200250	20X25 cm	5	24	
PADEX200300	PADEXC200300	PADEXALGI200300	20X30 cm	5	24	
PADEX200400	PADEXC200400	PADEXALGI200400	20X40 cm	5	24	
PADEX7575	PADEXC7575	PADEXALGI7575	7.5X7.5 cm	10	48	
PADEX200500	PADEXC200500	PADEXALGI200500	20x50 cm			
PADEX220220	PADEXC220220	PADEXALGI220220	22x22 cm			
PADEX250450	PADEXC250450	PADEXALGI250450	25x45 cm			
Cure -Aid <sup>®</sup> eXsorb <sup>®</sup>			Size	Pieces/ box	Boxes/ Case	Case Dimensions
CUREEX938			9x38 mm	100	60	
			9.5x38mm			
CUREEX1657			16x57mm	100	60	
CUREEX1938			19x38mm	100	60	
CUREEX1963			19x63mm	100	60	
CUREEX1964			19x64mm	100	60	
CUREEX1972			19x72mm	100	60	
CUREEX1976			19x76mm	100	60	
CUREEX2075			20x75mm	100	60	
CUREEX2572			25x72mm	100	60	
CUREEX2576			25x76mm	100	60	
CUREEX3055			30x55mm	100	60	
CUREEX3064			30x64mm	100	60	
CUREEX3072			30x72mm	100	60	
CUREEX3838			38x38mm	100	60	
CUREEXANCHOR			38x72mm Anchor	50	60	
CUREEX5172			51x72mm	50	60	
CUREEX50100			50x100mm	50	60	
CUREEX5060			50x60mm Fingertip	50	60	
CUREEX4563			45x63mm butterfly	50	60	
CUREEX2040			20x40 mm	100	60	
CUREEX2375			23x75mm	100	60	
CUREEX2525			25x25 mm	100	60	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

<b>Cure -Aid<sup>®</sup> eXsorb<sup>®</sup></b>	<b>Size</b>	<b>Pieces/ box</b>	<b>Boxes/ Case</b>	<b>Case Dimensions</b>
CUREEX2582	25x82 mm	100	60	
CUREEX2956	29x56 mm	100	60	
CUREEX3050	30x50 mm	100	60	
CUREEX4070	40x70 mm	50	60	
CUREEX5050	50x50 mm	100	60	
CUREEX5072	50x72 mm	50	60	
CUREEX2050	20x50 mm	100	60	
<b>Cure -Aid<sup>®</sup> Dressing Strip eXsorb<sup>®</sup></b>	<b>Size</b>	<b>Pieces/ box</b>	<b>Boxes/ Case</b>	<b>Case Dimensions</b>
DERMIEX2040	2x4 cm	100	60	
DERMIEX2060	2x6 cm	100	60	
DERMIEX20120	2X12cm	100	60	
DERMIEX30120	3x12cm	100	60	
DERMIEX20180	2x18cm	100	60	
DERMIEX30180	3x18cm	100	60	
DERMIEX40100	4X10cm	10	120	
DERMIEX60100	6X10cm	10	120	
DERMIEX70100	7x10cm	10	120	
DERMIEX80100	8x10cm	10	120	
DERMIEX601000	6x100cm	1	120	
DERMIEX801000	8x100cm	1	120	
<b>Pharmapore<sup>®</sup> eXsorb</b>	<b>Size</b>	<b>Pieces/ box</b>	<b>Boxes/ Case</b>	<b>Case Dimensions</b>
POREEX5070	5x7cm	10	48	
POREEX6070	6x7cm	10	48	
POREEX6080	6x8cm	10	48	
POREEX5090	5x9cm	10	48	
POREEX6083	6x8.3cm	10	48	
POREEX50100	5x10cm	10	48	
POREEX6082	6x8.25cm	10	48	
POREEX6085	6x8.5cm	10	48	
POREEX7080	7x8cm	10	48	
POREEX6580	6.5x8cm	10	48	
POREEX7575	7.5x7.5cm	10	48	
POREEX6090	6x9cm	10	48	
POREEX60100	6x10cm	10	48	
POREEX5082	5x8.25cm	10	48	
POREEX75100	7.5x10cm	10	48	
POREEX80100	8x10cm	10	48	
POREEX82100	8.25x10cm	10	48	
POREEX100100	10x10cm	10	48	
POREEX70120	7x12cm	10	48	
POREEX80120	8x12cm	10	48	
POREEX100120	10x12cm	10	48	
POREEX83120	8.3x12cm	10	48	
POREEX90150	9x15cm	10	48	
POREEX80150	8x15cm	10	48	
POREEX82150	8.25x15cm	10	48	
POREEX100140	10x14cm	10	48	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREEX100150	10x15cm	10	48	
POREEX150150	15x15cm	10	24	
POREEX83180	8.3x18cm	10	48	
POREEX82200	8.25x20cm	10	24	
POREEX90200	9x20cm	10	24	
POREEX100200	10x20cm	10	24	
POREEX200200	20x20cm	10	24	
POREEX83240	8.3x24cm	10	24	
POREEX82250	8.25x25cm	10	24	
POREEX90250	9x25cm	10	24	
POREEX100250	10x25cm	10	24	
POREEX82300	8.25x30cm	10	24	
POREEX90300	9x30cm	10	24	
POREEX100300	10x30cm	10	24	
POREEX90350	9x35cm	10	24	
POREEX82350	8.25x35cm	10	24	
POREEX100350	10x35cm	10	24	
POREEX83415	8.3x41.5cm	10	24	
POREEX100160	10x16 cm	10	48	
POREEX100210	10x21 cm	10	24	
POREEX100260	10x26 cm	10	24	
POREEX100400	10x40 cm	10	24	
POREEX120160	12x16 cm	10	48	
POREEX100380	10x38cm	10	24	
POREEX5072	5x7.2 cm	10	48	
POREEX5075	5x7.5 cm	10	48	
POREEX70100	7x10 cm	10	48	
POREEX6075	6x7.5 cm	10	48	
POREEX85155	8.5x15.5 cm	10	48	
POREEX85150	8.5x15cm	10	48	
POREEX8595	8.5x9.5 cm	10	48	
POREEX170170S	17X17 cm Sacral	10	24	
POREEX90100	9x10cm	10	48	
POREEX125125	12.5x12.5 cm			
POREEX150200	15x20 cm			
POREEX170170	17X17 cm			
Pharmapore <sup>®</sup> eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREEXC5070	5x7cm	10	48	
POREEXC6070	6x7cm	10	48	
POREEXC6080	6x8cm	10	48	
POREEXC5090	5x9cm	10	48	
POREEXC6083	6x8.3cm	10	48	
POREEXC50100	5x10cm	10	48	
POREEXC6082	6x8.25cm	10	48	
POREEXC6085	6x8.5cm	10	48	
POREEXC7080	7x8cm	10	48	
POREEXC6580	6.5x8cm	10	48	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREEXC7575	7.5x7.5cm	10	48	
POREEXC6090	6x9cm	10	48	
POREEXC60100	6x10cm	10	48	
POREEXC5082	5x8.25cm	10	48	
POREEXC75100	7.5x10cm	10	48	
POREEXC80100	8x10cm	10	48	
POREEXC82100	8.25x10cm	10	48	
POREEXC100100	10x10cm	10	48	
POREEXC70120	7x12cm	10	48	
POREEXC80120	8x12cm	10	48	
POREEXC100120	10x12cm	10	48	
POREEXC83120	8.3x12cm	10	48	
POREEXC90150	9x15cm	10	48	
POREEXC80150	8x15cm	10	48	
POREEXC82150	8.25x15cm	10	48	
POREEXC100140	10x14cm	10	48	
POREEXC100150	10x15cm	10	48	
POREEXC150150	15x15cm	10	24	
POREEXC83180	8.3x18cm	10	48	
POREEXC82200	8.25x20cm	10	24	
POREEXC90200	9x20cm	10	24	
POREEXC100200	10x20cm	10	24	
POREEXC200200	20x20cm	10	24	
POREEXC83240	8.3x24cm	10	24	
POREEXC82250	8.25x25cm	10	24	
POREEXC90250	9x25cm	10	24	
POREEXC100250	10x25cm	10	24	
POREEXC82300	8.25x30cm	10	24	
POREEXC90300	9x30cm	10	24	
POREEXC100300	10x30cm	10	24	
POREEXC90350	9x35cm	10	24	
POREEXC82350	8.25x35cm	10	24	
POREEXC100350	10x35cm	10	24	
POREEXC83415	8.3x41.5cm	10	24	
POREEXC100160	10x16 cm	10	48	
POREEXC100210	10x21 cm	10	24	
POREEXC100260	10x26 cm	10	24	
POREEXC100400	10x40 cm	10	24	
POREEXC120160	12x16 cm	10	48	
POREEXC100380	10x38cm	10	24	
POREEXC5072	5x7.2 cm	10	48	
POREEXC5075	5x7.5 cm	10	48	
POREEXC70100	7x10 cm	10	48	
POREEXC6075	6x7.5 cm	10	48	
POREEXC85155	8.5x15.5 cm	10	48	
POREEXC85150	8.5x15cm	10	48	
POREEXC5682	5.6x8.2cm	10	48	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREEXC8595	8.5x9.5 cm	10	48	
POREEXC170170S	17X17 CM Sacral	10	24	
POREEXC90100	9x10cm	10	48	
Pharmapore <sup>®</sup> PU eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUEX5070	5x7cm	10	48	
POREPUEX6070	6x7cm	10	48	
POREPUEX6080	6x8cm	10	48	
POREPUEX5090	5x9cm	10	48	
POREPUEX6083	6x8.3cm	10	48	
POREPUEX50100	5x10cm	10	48	
POREPUEX6082	6x8.25cm	10	48	
POREPUEX6085	6x8.5cm	10	48	
POREPUEX7080	7x8cm	10	48	
POREPUEX6580	6.5x8cm	10	48	
POREPUEX7575	7.5x7.5cm	10	48	
POREPUEX6090	6x9cm	10	48	
POREPUEX60100	6x10cm	10	48	
POREPUEX5082	5x8.25cm	10	48	
POREPUEX75100	7.5x10cm	10	48	
POREPUEX80100	8x10cm	10	48	
POREPUEX82100	8.25x10cm	10	48	
POREPUEX100100	10x10cm	10	48	
POREPUEX70120	7x12cm	10	48	
POREPUEX80120	8x12cm	10	48	
POREPUEX100120	10x12cm	10	48	
POREPUEX83120	8.3x12cm	10	48	
POREPUEX90150	9x15cm	10	48	
POREPUEX80150	8x15cm	10	48	
POREPUEX82150	8.25x15cm	10	48	
POREPUEX100140	10x14cm	10	48	
POREPUEX100150	10x15cm	10	48	
POREPUEX150150	15x15cm	10	24	
POREPUEX83180	8.3x18cm	10	48	
POREPUEX82200	8.25x20cm	10	24	
POREPUEX90200	9x20cm	10	24	
POREPUEX100200	10x20cm	10	24	
POREPUEX200200	20x20cm	10	24	
POREPUEX83240	8.3x24cm	10	24	
POREPUEX82250	8.25x25cm	10	24	
POREPUEX90250	9x25cm	10	24	
POREPUEX100250	10x25cm	10	24	
POREPUEX82300	8.25x30cm	10	24	
POREPUEX90300	9x30cm	10	24	
POREPUEX100300	10x30cm	10	24	
POREPUEX90350	9x35cm	10	24	
POREPUEX82350	8.25x35cm	10	24	
POREPUEX100350	10x35cm	10	24	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> PU eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUEx83415	8.3x41.5cm	10	24	
POREPUEx100160	10x16 cm	10	48	
POREPUEx100210	10x21 cm	10	24	
POREPUEx100260	10x26 cm	10	24	
POREPUEx100400	10x40 cm	10	24	
POREPUEx120160	12x16 cm	10	48	
POREPUEx100380	10x38cm	10	24	
POREPUEx5072	5x7.2 cm	10	48	
POREPUEx5075	5x7.5 cm	10	48	
POREPUEx70100	7x10 cm	10	48	
POREPUEx6075	6x7.5 cm	10	48	
POREPUEx85155	8.5x15.5 cm	10	48	
POREPUEx85150	8.5x15cm	10	48	
POREPUEx8595	8.5x9.5 cm	10	48	
POREPUEx170170S	17X17 cm Sacral	10	24	
POREPUEx90100	9x10cm	10	48	
POREPUEx125125	12.5x12.5 cm			
POREPUEx150200	15x20 cm			
POREPUEx170170	17X17 cm			
Pharmapore <sup>®</sup> PU eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUExC5070	5x7cm	10	48	
POREPUExC6070	6x7cm	10	48	
POREPUExC6080	6x8cm	10	48	
POREPUExC5090	5x9cm	10	48	
POREPUExC6083	6x8.3cm	10	48	
POREPUExC50100	5x10cm	10	48	
POREPUExC6082	6x8.25cm	10	48	
POREPUExC6085	6x8.5cm	10	48	
POREPUExC7080	7x8cm	10	48	
POREPUExC6580	6.5x8cm	10	48	
POREPUExC7575	7.5x7.5cm	10	48	
POREPUExC6090	6x9cm	10	48	
POREPUExC60100	6x10cm	10	48	
POREPUExC5082	5x8.25cm	10	48	
POREPUExC75100	7.5x10cm	10	48	
POREPUExC80100	8x10cm	10	48	
POREPUExC82100	8.25x10cm	10	48	
POREPUExC100100	10x10cm	10	48	
POREPUExC70120	7x12cm	10	48	
POREPUExC80120	8x12cm	10	48	
POREPUExC100120	10x12cm	10	48	
POREPUExC83120	8.3x12cm	10	48	
POREPUExC90150	9x15cm	10	48	
POREPUExC80150	8x15cm	10	48	
POREPUExC82150	8.25x15cm	10	48	
POREPUExC100140	10x14cm	10	48	
POREPUExC100150	10x15cm	10	48	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> PU eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUExc150150	15x15cm	10	24	
POREPUExc83180	8.3x18cm	10	48	
POREPUExc82200	8.25x20cm	10	24	
POREPUExc90200	9x20cm	10	24	
POREPUExc100200	10x20cm	10	24	
POREPUExc200200	20x20cm	10	24	
POREPUExc83240	8.3x24cm	10	24	
POREPUExc82250	8.25x25cm	10	24	
POREPUExc90250	9x25cm	10	24	
POREPUExc100250	10x25cm	10	24	
POREPUExc82300	8.25x30cm	10	24	
POREPUExc90300	9x30cm	10	24	
POREPUExc100300	10x30cm	10	24	
POREPUExc90350	9x35cm	10	24	
POREPUExc82350	8.25x35cm	10	24	
POREPUExc100350	10x35cm	10	24	
POREPUExc83415	8.3x41.5cm	10	24	
POREPUExc100160	10x16 cm	10	48	
POREPUExc100210	10x21 cm	10	24	
POREPUExc100260	10x26 cm	10	24	
POREPUExc100400	10x40 cm	10	24	
POREPUExc120160	12x16 cm	10	48	
POREPUExc100380	10x38cm	10	24	
POREPUExc5072	5x7.2 cm	10	48	
POREPUExc5075	5x7.5 cm	10	48	
POREPUExc70100	7x10 cm	10	48	
POREPUExc6075	6x7.5 cm	10	48	
POREPUExc85155	8.5x15.5 cm	10	48	
POREPUExc85150	8.5x15cm	10	48	
POREPUExc5682	5.6x8.2cm	10	48	
POREPUExc8595	8.5x9.5 cm	10	48	
POREPUExc170170S	17X17 CM Sacral	10	24	
POREPUExc90100	9x10cm	10	48	
Pharmapore <sup>®</sup> PU FS eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUFSEX5070	5x7cm	10	48	
POREPUFSEX6070	6x7cm	10	48	
POREPUFSEX6080	6x8cm	10	48	
POREPUFSEX5090	5x9cm	10	48	
POREPUFSEX6083	6x8.3cm	10	48	
POREPUFSEX50100	5x10cm	10	48	
POREPUFSEX6082	6x8.25cm	10	48	
POREPUFSEX6085	6x8.5cm	10	48	
POREPUFSEX7080	7x8cm	10	48	
POREPUFSEX6580	6.5x8cm	10	48	
POREPUFSEX7575	7.5x7.5cm	10	48	
POREPUFSEX6090	6x9cm	10	48	
POREPUFSEX60100	6x10cm	10	48	



# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> PU FS eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUFSEX5082	5x8.25cm	10	48	
POREPUFSEX75100	7.5x10cm	10	48	
POREPUFSEX80100	8x10cm	10	48	
POREPUFSEX82100	8.25x10cm	10	48	
POREPUFSEX100100	10x10cm	10	48	
POREPUFSEX70120	7x12cm	10	48	
POREPUFSEX80120	8x12cm	10	48	
POREPUFSEX100120	10x12cm	10	48	
POREPUFSEX83120	8.3x12cm	10	48	
POREPUFSEX90150	9x15cm	10	48	
POREPUFSEX80150	8x15cm	10	48	
POREPUFSEX82150	8.25x15cm	10	48	
POREPUFSEX100140	10x14cm	10	48	
POREPUFSEX100150	10x15cm	10	48	
POREPUFSEX150150	15x15cm	10	24	
POREPUFSEX83180	8.3x18cm	10	48	
POREPUFSEX82200	8.25x20cm	10	24	
POREPUFSEX90200	9x20cm	10	24	
POREPUFSEX100200	10x20cm	10	24	
POREPUFSEX200200	20x20cm	10	24	
POREPUFSEX83240	8.3x24cm	10	24	
POREPUFSEX82250	8.25x25cm	10	24	
POREPUFSEX90250	9x25cm	10	24	
POREPUFSEX100250	10x25cm	10	24	
POREPUFSEX82300	8.25x30cm	10	24	
POREPUFSEX90300	9x30cm	10	24	
POREPUFSEX100300	10x30cm	10	24	
POREPUFSEX90350	9x35cm	10	24	
POREPUFSEX82350	8.25x35cm	10	24	
POREPUFSEX100350	10x35cm	10	24	
POREPUFSEX83415	8.3x41.5cm	10	24	
POREPUFSEX100160	10x16 cm	10	48	
POREPUFSEX100210	10x21 cm	10	24	
POREPUFSEX100260	10x26 cm	10	24	
POREPUFSEX100400	10x40 cm	10	24	
POREPUFSEX120160	12x16 cm	10	48	
POREPUFSEX100380	10x38cm	10	24	
POREPUFSEX5072	5x7.2 cm	10	48	
POREPUFSEX5075	5x7.5 cm	10	48	
POREPUFSEX70100	7x10 cm	10	48	
POREPUFSEX6075	6x7.5 cm	10	48	
POREPUFSEX85155	8.5x15.5 cm	10	48	
POREPUFSEX85150	8.5x15cm	10	48	
POREPUFSEX8595	8.5x9.5 cm	10	48	
POREPUFSEX170170S	17X17 cm Sacral	10	24	
POREPUFSEX90100	9x10cm	10	48	
POREPUFSEX125125	12.5x12.5 cm			

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> PU FS eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUFSEX150200	15x20 cm			
POREPUFSEX170170	17X17 cm			
Pharmapore <sup>®</sup> PU FS eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREPUFSEXC5070	5x7cm	10	48	
POREPUFSEXC6070	6x7cm	10	48	
POREPUFSEXC6080	6x8cm	10	48	
POREPUFSEXC5090	5x9cm	10	48	
POREPUFSEXC6083	6x8.3cm	10	48	
POREPUFSEXC50100	5x10cm	10	48	
POREPUFSEXC6082	6x8.25cm	10	48	
POREPUFSEXC6085	6x8.5cm	10	48	
POREPUFSEXC7080	7x8cm	10	48	
POREPUFSEXC6580	6.5x8cm	10	48	
POREPUFSEXC7575	7.5x7.5cm	10	48	
POREPUFSEXC6090	6x9cm	10	48	
POREPUFSEXC60100	6x10cm	10	48	
POREPUFSEXC5082	5x8.25cm	10	48	
POREPUFSEXC75100	7.5x10cm	10	48	
POREPUFSEXC80100	8x10cm	10	48	
POREPUFSEXC82100	8.25x10cm	10	48	
POREPUFSEXC100100	10x10cm	10	48	
POREPUFSEXC70120	7x12cm	10	48	
POREPUFSEXC80120	8x12cm	10	48	
POREPUFSEXC100120	10x12cm	10	48	
POREPUFSEXC83120	8.3x12cm	10	48	
POREPUFSEXC90150	9x15cm	10	48	
POREPUFSEXC80150	8x15cm	10	48	
POREPUFSEXC82150	8.25x15cm	10	48	
POREPUFSEXC100140	10x14cm	10	48	
POREPUFSEXC100150	10x15cm	10	48	
POREPUFSEXC150150	15x15cm	10	24	
POREPUFSEXC83180	8.3x18cm	10	48	
POREPUFSEXC82200	8.25x20cm	10	24	
POREPUFSEXC90200	9x20cm	10	24	
POREPUFSEXC100200	10x20cm	10	24	
POREPUFSEXC200200	20x20cm	10	24	
POREPUFSEXC83240	8.3x24cm	10	24	
POREPUFSEXC82250	8.25x25cm	10	24	
POREPUFSEXC90250	9x25cm	10	24	
POREPUFSEXC100250	10x25cm	10	24	
POREPUFSEXC82300	8.25x30cm	10	24	
POREPUFSEXC90300	9x30cm	10	24	
POREPUFSEXC100300	10x30cm	10	24	
POREPUFSEXC90350	9x35cm	10	24	
POREPUFSEXC82350	8.25x35cm	10	24	
POREPUFSEXC100350	10x35cm	10	24	
POREPUFSEXC83415	8.3x41.5cm	10	24	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> PU FS eXsorb Carbon		Size	Pieces/ box	Boxes/ Case	Case Dimensions	
	POREPUFSEXC100160	10x16 cm	10	48		
	POREPUFSEXC100210	10x21 cm	10	24		
	POREPUFSEXC100260	10x26 cm	10	24		
	POREPUFSEXC100400	10x40 cm	10	24		
	POREPUFSEXC120160	12x16 cm	10	48		
	POREPUFSEXC100380	10x38cm	10	24		
	POREPUFSEXC5072	5x7.2 cm	10	48		
	POREPUFSEXC5075	5x7.5 cm	10	48		
	POREPUFSEXC70100	7x10 cm	10	48		
	POREPUFSEXC6075	6x7.5 cm	10	48		
	POREPUFSEXC85155	8.5x15.5 cm	10	48		
	POREPUFSEXC85150	8.5x15cm	10	48		
	POREPUFSEXC8595	8.5x9.5 cm	10	48		
	POREPUFSEXC90100	9x10cm	10	48		
	POREPUFSEXC170170S	17X17 CM Sacral	10	24		
Pharmapore <sup>®</sup> IV eXsorb	Pharmapore <sup>®</sup> PU I.V eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions	
	IVNWEX6070	IVPUEX6070	6X7CM	100	16	
	IVNWEX7070	IVPUEX7070	7X7CM	100	16	
	IVNWEX7080	IVPUEX7080	7X8CM	100	16	
	IVNWEX80100	IVPUEX80100	8X10CM	100	16	
	IVNWEX82100	IVPUEX82100	8.25X10CM	100	16	
	IVNWEX6085	IVPUEX6085	6X8.5CM	100	16	
	IVNWEX7085	IVPUEX7085	7X8.5CM	100	16	
	IVNWEX85100	IVPUEX85100	8.5X10CM	100	16	
	IVNWEX7090	IVPUEX7090	7X9CM	100	16	
	IVNWEX90100	IVPUEX90100	9X10CM	100	16	
	IVNWEX100120	IVPUEX100120	10X12CM	60	16	
	IVNWEX100140	IVPUEX100140	10X14CM	60	16	
	IVNWEX6080	IVPUEX6080	6X8CM	100	16	
	IVNWEX85115	IVPUEX85115	8.5X11.5CM	100	16	
	IVNWEX5057	IVPUEX5057	5X5.7CM	100	16	
	IVNWEX100155	IVPUEX100155	10x15.5 cm	60	16	
	IVNWEX5055	IVPUEX5055	5X5.5 cm	100	16	
	IVNWEX5060	IVPUEX5060	5X6 cm	100	16	
	IVNWEX6570	IVPUEX6570	6.5X7cm	100	16	
Pharmapore <sup>®</sup> Ultra eXsorb		Size	Pieces/ box	Boxes/ Case	Case Dimensions	
	POREUEX5070	5x7cm	10	48		
	POREUEX6070	6x7cm	10	48		
	POREUEX6080	6x8cm	10	48		
	POREUEX5090	5x9cm	10	48		
	POREUEX6083	6x8.3cm	10	48		
	POREUEX50100	5x10cm	10	48		
	POREUEX6082	6x8.25cm	10	48		
	POREUEX6085	6x8.5cm	10	48		
	POREUEX7080	7x8cm	10	48		
	POREUEX6580	6.5x8cm	10	48		
	POREUEX7575	7.5x7.5cm	10	48		

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> Ultra eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREUEX6090	6x9cm	10	48	
POREUEX60100	6x10cm	10	48	
POREUEX5082	5x8.25cm	10	48	
POREUEX75100	7.5x10cm	10	48	
POREUEX80100	8x10cm	10	48	
POREUEX82100	8.25x10cm	10	48	
POREUEX100100	10x10cm	10	48	
POREUEX70120	7x12cm	10	48	
POREUEX80120	8x12cm	10	48	
POREUEX100120	10x12cm	10	48	
POREUEX83120	8.3x12cm	10	48	
POREUEX90150	9x15cm	10	48	
POREUEX80150	8x15cm	10	48	
POREUEX82150	8.25x15cm	10	48	
POREUEX100140	10x14cm	10	48	
POREUEX100150	10x15cm	10	48	
POREUEX150150	15x15cm	10	24	
POREUEX83180	8.3x18cm	10	48	
POREUEX82200	8.25x20cm	10	24	
POREUEX90200	9x20cm	10	24	
POREUEX100200	10x20cm	10	24	
POREUEX200200	20x20cm	10	24	
POREUEX83240	8.3x24cm	10	24	
POREUEX82250	8.25x25cm	10	24	
POREUEX90250	9x25cm	10	24	
POREUEX100250	10x25cm	10	24	
POREUEX82300	8.25x30cm	10	24	
POREUEX90300	9x30cm	10	24	
POREUEX100300	10x30cm	10	24	
POREUEX90350	9x35cm	10	24	
POREUEX82350	8.25x35cm	10	24	
POREUEX100350	10x35cm	10	24	
POREUEX83415	8.3x41.5cm	10	24	
POREUEX100160	10x16 cm	10	48	
POREUEX100210	10x21 cm	10	24	
POREUEX100260	10x26 cm	10	24	
POREUEX100400	10x40 cm	10	24	
POREUEX120160	12x16 cm	10	48	
POREUEX100380	10x38cm	10	24	
POREUEX5072	5x7.2 cm	10	48	
POREUEX5075	5x7.5 cm	10	48	
POREUEX70100	7x10 cm	10	48	
POREUEX6075	6x7.5 cm	10	48	
POREUEX85155	8.5x15.5 cm	10	48	
POREUEX85150	8.5x15cm	10	48	
POREUEX8595	8.5x9.5 cm	10	48	
POREUEX170170S	17X17 cm Sacral	10	24	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family


Pharmapore <sup>®</sup> Ultra eXsorb	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREUEX90100	9x10cm	10	48	
POREUEX125125	12.5x12.5 cm			
POREUEX150200	15x20 cm			
POREUEX170170	17X17 cm			
Pharmapore <sup>®</sup> Ultra eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREUEXC5070	5x7cm	10	48	
POREUEXC6070	6x7cm	10	48	
POREUEXC6080	6x8cm	10	48	
POREUEXC5090	5x9cm	10	48	
POREUEXC6083	6x8.3cm	10	48	
POREUEXC50100	5x10cm	10	48	
POREUEXC6082	6x8.25cm	10	48	
POREUEXC6085	6x8.5cm	10	48	
POREUEXC7080	7x8cm	10	48	
POREUEXC6580	6.5x8cm	10	48	
POREUEXC7575	7.5x7.5cm	10	48	
POREUEXC6090	6x9cm	10	48	
POREUEXC60100	6x10cm	10	48	
POREUEXC5082	5x8.25cm	10	48	
POREUEXC75100	7.5x10cm	10	48	
POREUEXC80100	8x10cm	10	48	
POREUEXC82100	8.25x10cm	10	48	
POREUEXC100100	10x10cm	10	48	
POREUEXC70120	7x12cm	10	48	
POREUEXC80120	8x12cm	10	48	
POREUEXC100120	10x12cm	10	48	
POREUEXC83120	8.3x12cm	10	48	
POREUEXC90150	9x15cm	10	48	
POREUEXC80150	8x15cm	10	48	
POREUEXC82150	8.25x15cm	10	48	
POREUEXC100140	10x14cm	10	48	
POREUEXC100150	10x15cm	10	48	
POREUEXC150150	15x15cm	10	24	
POREUEXC83180	8.3x18cm	10	48	
POREUEXC82200	8.25x20cm	10	24	
POREUEXC90200	9x20cm	10	24	
POREUEXC100200	10x20cm	10	24	
POREUEXC200200	20x20cm	10	24	
POREUEXC83240	8.3x24cm	10	24	
POREUEXC82250	8.25x25cm	10	24	
POREUEXC90250	9x25cm	10	24	
POREUEXC100250	10x25cm	10	24	
POREUEXC82300	8.25x30cm	10	24	
POREUEXC90300	9x30cm	10	24	
POREUEXC100300	10x30cm	10	24	
POREUEXC90350	9x35cm	10	24	
POREUEXC82350	8.25x35cm	10	24	

# eXsorb<sup>®</sup> and eXsorb<sup>®</sup> Carbon Family

Pharmapore <sup>®</sup> Ultra eXsorb Carbon	Size	Pieces/ box	Boxes/ Case	Case Dimensions
POREUEXC100350	10x35cm	10	24	
POREUEXC83415	8.3x41.5cm	10	24	
POREUEXC100160	10x16 cm	10	48	
POREUEXC100210	10x21 cm	10	24	
POREUEXC100260	10x26 cm	10	24	
POREUEXC100400	10x40 cm	10	24	
POREUEXC120160	12x16 cm	10	48	
POREUEXC100380	10x38cm	10	24	
POREUEXC5072	5x7.2 cm	10	48	
POREUEXC5075	5x7.5 cm	10	48	
POREUEXC70100	7x10 cm	10	48	
POREUEXC6075	6x7.5 cm	10	48	
POREUEXC85155	8.5x15.5 cm	10	48	
POREUEXC85150	8.5x15cm	10	48	
POREUEXC5682	5.6x8.2cm	10	48	
POREUEXC8595	8.5x9.5 cm	10	48	
POREUEXC170170S	17X17 CM Sacral	10	24	
POREUEXC90100	9x10cm	10	48	

- Amria - Alexandria - Egypt  
Tel. : 002 03 4500264 Fax : 002 03 4500263  
- Kafr-El-Zayat P.O.Box 38 Egypt  
Tel. : 002 040 2545643 Fax : 002 040 2545645  
E-mail : sales@pharmaplast-online.com

 **Pharmaplast**  
www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3 )						
Customer : Pharmaplast		Product name : eXsorb and eXsorb Carbon Family						
Issue date : 16-11-2016		Prepared by : Marina Medhat						
Colors : 		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	04-02-2019	Essential requirement	changes in able of size by haidy samir					
3	04/04/2019	-Text changed upon request from ( Haidy Samir)	To be comply with ISO 15223-1. Due to the brexit issues.					
4	05/05/2019	- Text changed upon request from ( Jostina Ramsis)	To be comply with ISO 15223-1. Due to the brexit issues.					
5	10-01-2022	- Updating table of sizes	File Maintenance					