

ESPUMA® Gentle Family

Soft Silicone Polyurethane Foam Dressings



Description:

ESPUMA® Gentle is absorbent wound dressing made of polyurethane foam .The outermost layer is bacteria-proof, virus-proof and water-proof polyurethane film. The wound contact layer is perforated soft silicone gel that allows painless and atraumatic removal and permits the passage of exudate to the absorbent foam layer.

ESPUMA® Gentle Lite consists of thin absorbent polyurethane foam layer.

ESPUMA® Gentle Comfort consists of absorbent foam layer having an adhesive outermost polyurethane border, so the dressing can be applied without fixation tape or secondary dressing.

ESPUMA® Gentle Lite Comfort has a thin foam layer with an adhesive outermost PU border.

ESPUMA® Gentle C Comfort is an absorbent foam with additional non-woven layer containing activated carbon to adsorb offensive odour and with an adhesive outermost PU border.

ESPUMA® Gentle Lite C Comfort differs from the previous variant in having a thinner polyurethane foam layer.

ESPUMA® Gentle Xtra Comfort is a multi-layered dressing consisting of absorbent foam with additional highly absorbent layer and having an adhesive outermost border.

ESPUMA® Gentle Lite Xtra Comfort differs from the previous variant in having a thinner polyurethane foam layer.

ESPUMA® Gentle C Xtra Comfort is a multi-layered dressing consisting of absorbent foam, with additional highly absorbent layer plus a non-woven layer containing activated carbon to adsorb offensive odour and having an adhesive outermost border.

ESPUMA® Gentle Lite C Xtra Comfort differs from the previous variant in having a thinner polyurethane foam layer.

ESPUMA® Gentle Transfer is composed of thin polyurethane foam layer, but without polyurethane top sheet, so that the dressing can transfer exudates away from the wound into a secondary absorbent dressing.

Intended use:

ESPUMA® Gentle Family is an effective tool for moist wound healing, having high capacity of autolytic debridement and it is particularly useful in preventing dressing-related trauma especially in people with fragile skin. It's indicated for a wide range of exuding wounds as: pressure ulcers, leg ulcers, diabetic foot ulcers, skin graft, first and second degree burns. Can be used in combination with compression therapy.

The Lite variants are indicated for low exuding and traumatic wounds as: abrasions, blisters and skin tears.

The variants containing carbon are for malodorous acute and chronic wounds as carbon adsorb offensive odor.

Features:

- Comfortable: as the dressing is easy to remove and minimize trauma during dressing changes.
- Cushioning effect: as it acts as a barrier against irritation and constant pressure on the skin.
- Breathable: the semi-permeable backing allows moisture to escape.
- Prevents skin maceration: as the perforated silicone film allows free passage of wound exudate to the foam.
- Additional carbon layer in some variants: adsorbs offensive odour in case of malodorous wounds.
- Additional superabsorbent layer in some variants: increases the absorption capacity of the dressing for highly exuding wounds.

ESPUMA® Gentle Family

Soft Silicone Polyurethane Foam Dressings

| ESPUMA® Gentle | ESPUMA® Gentle Lite | ESPUMA® Gentle Transfer | Size | Units/Box | Boxes / Case | Case dimensions |
|--|---|---|----------------|-----------|--------------|-----------------|
| ESPG7575 | ESPGL7575 | FoamgentleT7575 | □ 75 x 75 mm | 10 | 48 | |
| ESPG75100 | ESPGL75100 | FoamgentleT75100 | □ 75 x 100 mm | 10 | 48 | |
| ESPG100100 | ESPGL100100 | FoamgentleT100100 | □ 100 x 100 mm | 5 | 24 | 47x47x41 cm |
| ESPG100150 | ESPGL100150 | FoamgentleT100150 | □ 150 x 150 mm | 5 | 24 | 43x38x25 cm |
| ESPG100200 | ESPGL100200 | FoamgentleT100200 | □ 100 x 200 mm | 5 | 24 | |
| ESPG100250 | ESPGL100250 | FoamgentleT100250 | □ 100 x 250 mm | 5 | 24 | |
| ESPG100300 | ESPGL100300 | FoamgentleT100300 | □ 100 x 300 mm | 5 | 24 | |
| ESPG125125 | ESPGL125125 | FoamgentleT125125 | □ 125 x 125 mm | 5 | 24 | |
| ESPG130200 | ESPGL130200 | FoamgentleT130200 | ∞ 130 x 200 mm | 5 | 24 | |
| ESPG150150 | ESPGL150150 | FoamgentleT150150 | □ 150 x 150 mm | 5 | 24 | |
| ESPG175175 | ESPGL175175 | FoamgentleT175175 | □ 175 x 175 mm | 3 | 24 | |
| ESPG150200 | ESPGL150200 | FoamgentleT150200 | □ 150 x 200 mm | 3 | 24 | 45x26x30 cm |
| ESPG200200 | ESPGL200200 | FoamgentleT200200 | □ 200 x 200 mm | 3 | 24 | 60x28x32 cm |
| ESPG200300 | ESPGL200300 | FoamgentleT200300 | □ 200 x 300 mm | 3 | 24 | |
| ESPG200500 | ESPGL200500 | FoamgentleT200500 | □ 200 x 500 mm | 3 | 24 | |
| ESPG 6080 <small>Tracheostomy</small> | ESPGL6080 <small>Tracheostomy</small> | FoamgentleT <small>Tracheostomy</small> | ∞ 60 x 80 mm | 10 | 48 | 58x47x41 cm |
| ESPG6565 <small>Tracheostomy</small> | ESPGL6565 <small>Tracheostomy</small> | FoamgentleT <small>Tracheostomy</small> | ∞ 65 x 65 mm | 10 | 48 | |
| ESPG8080 <small>Tracheostomy</small> | ESPGL8080 <small>Tracheostomy</small> | FoamgentleT <small>Tracheostomy</small> | ∞ 80 x 80 mm | 10 | 48 | |
| ESPG8090 <small>Tracheostomy</small> | ESPGL8090 <small>Tracheostomy</small> | FoamgentleT <small>Tracheostomy</small> | ∞ 80 x 90 mm | 10 | 48 | |
| ESPG165180 <small>Concave</small> | ESPGL165180 <small>Concave</small> | FoamgentleT <small>Concave</small> | ◐ 165 x 180 mm | 3 | 24 | |
| ESPG170170 <small>Small Sacral</small> | ESPGL170170 <small>Small Sacral</small> | FoamgentleT <small>Small Sacral</small> | ◑ 170 x 170 mm | 3 | 24 | |
| ESPG220220 <small>Sacral</small> | ESPGL 220220 <small>Sacral</small> | FoamgentleT <small>Sacral</small> | ◒ 220 x 220 mm | 3 | 24 | |

| ESPUMA® Gentle Comfort | ESPUMA® Gentle Lite Comfort | Size | Units/Box | Boxes / Case | Case dimensions |
|--|---|--------------|-----------|--------------|-----------------|
| ESPGC7575 | ESPGCL 7575 | 75 x 75 mm | 10 | 48 | |
| ESPGC75100 | ESPGCL 75100 | 75 x 100 mm | 10 | 48 | |
| ESPGC100100 | ESPGCL 100100 | 100 x 100 mm | 10 | 48 | 47x47x41 cm |
| ESPGC100150 | ESPGCL 100150 | 100 x 150 mm | 5 | 24 | |
| ESPGC100200 | ESPGCL 100200 | 100 x 200 mm | 5 | 24 | |
| ESPGC100250 | ESPGCL 100250 | 100 x 250 mm | 5 | 24 | |
| ESPGC100300 | ESPGCL 100300 | 100 x 300 mm | 5 | 24 | |
| ESPGC125125 | ESPGCL125125 | 125 x 125 mm | 5 | 24 | |
| ESPGC 150150 | ESPGCL 150150 | 150 x 150 mm | 5 | 24 | 43x38x25 cm |
| ESPGC 150200 | ESPGCL 150200 | 150 x 200 mm | 3 | 24 | 45x26x30 cm |
| ESPGC200200 | ESPGCL 200200 | 200 x 200 mm | 3 | 24 | 60x28x32 cm |
| ESPGC200300 | ESPGCL 200300 | 200 x 300 mm | 3 | 24 | |
| ESPGC200500 | ESPGCL 200500 | 200 x 500 mm | 3 | 24 | |
| ESPGC 165180 <small>Concave</small> | ESPGCL 165180 <small>Concave</small> | 165 x 180 mm | 5 | 24 | |
| ESPGC 170170 <small>Small Sacral</small> | ESPGCL 170170 <small>Small Sacral</small> | 170 x 170 mm | 3 | 24 | |
| ESPGC220220 <small>Sacral</small> | ESPGCL 220220 <small>Sacral</small> | 220 x 220 mm | 3 | 24 | |

| ESPUMA® Gentle C Comfort | ESPUMA® Gentle Lite C Comfort | Size | Units/Box | Boxes / Case | Case dimensions |
|---|--|--------------|-----------|--------------|-----------------|
| ESPGCC 7575 | ESPGCCL 7575 | 75 x 75 mm | 10 | 48 | 46x46x39 cm |
| ESPGCC 75100 | ESPGCCL 75100 | 75 x 100 mm | 10 | 48 | 46x46x39 cm |
| ESPGCC 100100 | ESPGCCL 100100 | 100 x 100 mm | 10 | 48 | 47x47x41 cm |
| ESPGCC 100150 | ESPGCCL 100150 | 100 x 150 mm | 5 | 24 | |
| ESPGCC 100200 | ESPGCCL 100200 | 100 x 200 mm | 5 | 24 | |
| ESPGCC 100250 | ESPGCCL 100250 | 100 x 250 mm | 5 | 24 | |
| ESPGCC 100300 | ESPGCCL 100300 | 100 x 300 mm | 5 | 24 | |
| ESPGCC125125 | ESPGCCL125125 | 125 x 125 mm | 5 | 24 | |
| ESPGCC 150150 | ESPGCCL 150150 | 150 x 150 mm | 5 | 24 | 43x38x25 cm |
| ESPGCC 150200 | ESPGCCL 150200 | 150 x 200 mm | 3 | 24 | 45x26x30 cm |
| ESPGCC200200 | ESPGCCL 200200 | 200 x 200 mm | 3 | 24 | 53x26x30 cm |
| ESPGCC200300 | ESPGCCL 200300 | 200 x 300 mm | 3 | 24 | |
| ESPGCC200500 | ESPGCCL 200500 | 200 x 500 mm | 3 | 24 | |
| ESPGCC 165180 <small>Concave</small> | ESPGCCL 165180 <small>Concave</small> | 165 x 180 mm | 5 | 24 | 54x26x33 cm |
| ESPGCC 170170 <small>Small Sacral</small> | ESPGCCL 170170 <small>Small Sacral</small> | 170 x 170 mm | 3 | 24 | |
| ESPGCC220220 <small>Sacral</small> | ESPGCCL 220220 <small>Sacral</small> | 220 x 220 mm | 3 | 24 | |

ESPUMA® Gentle Family

Soft Silicone Polyurethane Foam Dressings

| ESPUMA® Gentle Xtra Comfort | ESPUMA® Gentle Lite Xtra Comfort | Size | Units/Box | Boxes / Case | Case dimensions |
|--|---|--------------|-----------|--------------|-----------------|
| ESPGCEX 7575 | ESPGCEXL 7575 | 75 x 75 mm | 10 | 48 | 46x46x39 cm |
| ESPGCEX 75100 | ESPGCEXL 75100 | 75 x 100 mm | 10 | 48 | 46x46x39 cm |
| ESPGCEX 100100 | ESPGCEXL 100100 | 100 x 100 mm | 10 | 48 | 47x47x41 cm |
| ESPGCEX 100150 | ESPGCEXL 100150 | 100 x 150 mm | 5 | 24 | |
| ESPGCEX 100200 | ESPGCEXL 100200 | 100 x 200 mm | 5 | 24 | |
| ESPGCEX 100250 | ESPGCEXL 100250 | 100 x 250 mm | 5 | 24 | |
| ESPGCEX 100300 | ESPGCEXL 100300 | 100 x 300 mm | 5 | 24 | |
| ESPGCEX125125 | ESPGCEXL125125 | 125 x 125 mm | 5 | 24 | |
| ESPGCEX 150150 | ESPGCEXL 150150 | 150 x 150 mm | 5 | 24 | 43x38x25 cm |
| ESPGCEX 150200 | ESPGCEXL 150200 | 150 x 200 mm | 3 | 24 | 45x26x30 cm |
| ESPGCEX 200200 | ESPGCEXL 200200 | 200 x 200 mm | 3 | 24 | 53x26x30 cm |
| ESPGCEX 200300 | ESPGCEXL 200300 | 200 x 300 mm | 3 | 24 | |
| ESPGCEX 200500 | ESPGCEXL 200500 | 200 x 500 mm | 3 | 24 | |
| ESPGCEX 165180 <small>Concave</small> | ESPGCEXL 165180 <small>Concave</small> | 165 x 180 mm | 5 | 24 | 54x26x33 cm |
| ESPGCEX 170170 <small>Small Sacral</small> | ESPGCEXL 170170 <small>Small Sacral</small> | 170 x 170 mm | 3 | 24 | |
| ESPGCEX 220220 <small>Sacral</small> | ESPGCEXL 220220 <small>Sacral</small> | 220 x 220 mm | 3 | 24 | |

Other sizes available upon request

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
 www.pharmaplast-online.com

| Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com | | Brochure artwork template no: I-TM /AR/001-6 | | | | |
|---|---------------|--|----------------------|-----------------|----------------------|---------------|
| Customer : Pharmaplast | | Product name : Espuma Gentle | | | | |
| Issue date : 02-10-2012 | | Prepared by : Nader Nabil Zaki | | | | |
| Colors : <input type="text"/> | | Dimension : A4 | | | | |
| Issue Number | Revision Date | Description of Change | Rationale For Change | Approved by | | |
| | | | | Artwork manager | Marketing department | RA Department |
| 1 | N/A | New Issue | N/A | | | |
| 2 | 05-05-2019 | New Text by Jostina Ramsis | | | | |