

# Kemagel® Conductive



## Description:

Non-sterile water based non-irritating hydrogel.

## Mode of Action:

Kemagel® Conductive is a coupling media, in the form of water soluble gel, prevents reflection of the waves away at the soft tissue/air interface, by removing air from between the patient and transducer. It also reduces static when using a Doppler, especially for fetal development images.

## Features:

1. Water soluble and non-greasy, so it can be wiped off easily.
2. Its unique cohesive properties allow it to spread over the skin evenly and easily.
3. Will not stain clothing or damage transducers.
4. No fragrance.
5. No formaldehydes.
6. Non-irritant formula.
7. Cost-effective and easy to use.
8. Safe for both patients and machinery.
9. High diagnostic sensitivity.

## Intended Purpose:

For external use only. **Kemagel® Conductive** is intended to be used on intact skin where ultrasound scan is needed.

Product	Code	Size	PCS/case
<b>Kemagel® Conductive</b>	<b>GELCON6</b>	<b>6gm Sachet</b>	<b>2000</b>
<b>Kemagel® Conductive</b>	<b>GELCON250</b>	<b>250ml</b>	<b>48</b>
<b>Kemagel® Conductive</b>	<b>GELCON20</b>	<b>20 gm Sachet</b>	<b>1000</b>
<b>Kemagel® Conductive</b>	<b>GELCON1000</b>	<b>1000ml</b>	<b>6</b>
<b>Kemagel® Conductive</b>	<b>GELCON4000</b>	<b>4000ml</b>	<b>1</b>
<b>Kemagel® Conductive</b>	<b>GELCON5000</b>	<b>5000ml</b>	<b>2</b>
<b>Kemagel® Conductive</b>	<b>GELCON10000</b>	<b>10000ml</b>	<b>1</b>
<b>Kemagel® Conductive</b>	<b>GELCON15000</b>	<b>15000ml</b>	<b>1</b>
<b>Kemagel® Conductive</b>	<b>GELCON20000</b>	<b>20000ml</b>	<b>1</b>

- Amria - Alexandria - Egypt  
Tel. : 002 03 4500264 Fax : 002 03 4500263  
- Kafr-El-Zayat P.O.Box 38 Egypt  
Tel. : 002 040 2545643 Fax : 002 040 2545645  
E-mail : sales@pharmaplast-online.com

 **Pharmaplast**  
www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3 )						
Customer : Pharmaplast		Product name : Kemagel Conductive						
Issue date : 15-06-2016		Prepared by : Marina Medhat						
Colors : CMYK		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	02-02-2020	Adding new sizes to table of sizes	Marketing point of view					
3	28-04-2020	Amending in table of size according to CRS	Marketing point of view					
4	06-05-2020	New text by Rania Raaft	file maintenance					
5	18-05-2020	changing packaging configuration according to CRS	to reflect current packaging					
6	05-06-2022	Updating Codes	File Maintenance					
7	28-03-2026	adding new size 1000 ml	File Maintenance	Nader Nabil zaki		Ramy Edward	Merna Samuel	Aksany Gasser