


Pharmastrip® and Pharmastrip® Fabric



Description:

Pharmastrip®

Flexible skin closure strips made of a soft non-woven backing laminated with bacterial proof polyurethane layer coated with a pressure-sensitive, hypoallergenic adhesive. Can be used on contoured body areas or where bending and movement occur.

Pharmastrip® Fabric

Flexible skin closure strips made of non-elastic backing coated with a Pressure-sensitive, hypoallergenic adhesive. Used in conjunction with skin suture and staples or after their removal for wound support.

Features:

Pharmastrip®:

- **Provides a perfect wound support:** due to the unique structure of nonwoven with PU layer which gives flexibility and strength.
- **Breathable, waterproof and bacteria proof:** decreasing the risk of cross infection.
- **Prevents roll-up:** due to the smooth rounded edges.
- Tan color for a cosmetically pleasing appearance.
- Fast and easy application to save time.
- Available in sealed pouches E.O. sterilized.

Pharmastrip® Fabric:

- **Eliminates puncture marks:** for better cosmetic result than sutures or staples.
- **Breathable:** Permeable to moisture and air which promotes healing.
- Fast and easy application to save time.
- Available in sealed pouches E.O. sterilized.

Intended Purpose:

- Adhesive skin closure strips are used in the treatment of lacerations and surgical incisions.
- Can be used in combination with skin sutures and staples or after their removal for wound support.


Pharmastrip® and Pharmastrip® Fabric

Product Code	Pharmastrip®	Size	Pouches/box	Boxes/case	Case dimensions
PHST638	Pouch of 6 strips	6mm×38mm	200	16	56x36x25
PHST375	Pouch of 10 strips	3mm×75mm	100	16	56x36x25
PHST675	Pouch of 6 strips	6mm×75mm	100	16	56x36x25
PHST6100	Pouch of 10 strips	6mm×100mm	100	16	56x36x25
PHST12100	Pouch of 6 strips	12mm×100mm	100	16	56x36x25
PHST25127	Pouch of 6 strips	25mm×127mm	70	16	56x36x25

Product Code	Pharmastrip® Fabric	Size	Pouches/box	Boxes/case	Case dimensions
PHST438	Pouch of 8 strips	4mm×38mm	200	16	47x47x21
PHST476	Pouch of 8 strips	4mm×76mm	100	16	47x47x21
PHST6476	Pouch of 6 strips	6.4mm×76mm	100	16	47x47x21
PHST64102	Pouch of 10 strips	6.4mm×102mm	100	16	47x47x21
PHST13102	Pouch of 6 strips	13mm×102mm	100	16	47x47x21
PHST26102	Pouch of 4 strips	26mm×102mm	100	16	47x47x21

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast			Product name : pharmastrip and pharmastrip fabric					
Issue date : 02-05-2019			Prepared by : Marina Medhat					
Colors : 			Dimension : A4					
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	30-03-2020	adding case dimension to table of sizes	file maintenance					
3	09-01-2022	Updating table of size to be aligned with the new coding system	file maintenance					
4	07-03-2026	Updating text	file maintenance	Nada Youssef				